

FORA COVID-19 Antigen Rapid Test

ver 1.0 202011
312-4531100-005

INTENDED USE

The FORA COVID-19 Antigen Rapid Test is a lateral flow chromatographic immunoassay intended for the qualitative detection of nucleocapsid protein antigen from SARS-CoV-2 in fresh nasopharyngeal swab specimens directly from individuals who are suspected of COVID-19 by their healthcare provider.

The FORA COVID-19 Antigen Rapid Test is intended for use by trained clinical laboratory personnel specifically instructed and trained in *in vitro* diagnostic procedures and individuals trained in point of care settings.

The FORA COVID-19 Antigen Rapid Test does not differentiate between SARS-CoV and SARS-CoV-2.

Results are for the identification of the SARS-CoV-2 nucleocapsid protein antigen. The antigen is generally detectable in upper respiratory specimens during the acute phase of infection. Positive results indicate the presence of viral antigens, but clinical correlation with patient history and other diagnostic information is necessary to determine the status of the infection. Positive results do not rule out bacterial infection or co-infection with other viruses. The agent detected may not be the definite cause of the disease.

Negative results do not rule out SARS-CoV-2 infection and should not be used as the sole basis for treatment or patient management decisions, including infection control decisions. Negative results should be considered in the context of a patient's recent exposures, history and the presence of clinical signs and symptoms consistent with COVID-19, and confirmed with a molecular assay if necessary for patient management.

SUMMARY AND EXPLANATION

The coronavirus disease 2019 (COVID-19) is caused by the severe acute respiratory syndrome coronavirus 2 (SARS-CoV-2). SARS-CoV-2 is a coronavirus identified as the cause of an outbreak of respiratory illnesses first detected in Wuhan, China. The WHO declared on March 11, 2020 that COVID-19 was a pandemic.^[1] COVID-19 has caused millions of confirmed cases worldwide, including hundreds of thousands of deaths, and statistics are increasing.^[2] It has been reported that symptoms ranging from mild to severe may appear 2-14 days after exposure to SARS-CoV-2. People with these symptoms may suffer from COVID-19: fever, cough, and shortness of breath.^[3]

There are three main categories of COVID-19 tests: (1) The antibody test results (such as IgM and IgG) can show whether a person has been infected before. The human immune system may take 15 days in most patients to produce antibodies after infection.^[4] (2) The antigen test results can show whether there is currently a virus infection. According to CLSI standard of viral culture, the median incubation time is estimated to be 5.1 days with symptoms expected to be present within 12 days of infection.^[5] FORA COVID-19 Antigen Rapid Test uses the monoclonal antibody that specifically binds to the nucleocapsid (N) protein to determine the presence of the SARS-CoV-2 antigen. (3) Real-Time RT-PCR (Reverse-Transcription Polymerase Chain Reaction) is a molecular method intended for the qualitative detection of nucleic acid from SARS-CoV-2. SARS-CoV-2 RNA is converted into complementary DNA by using a reverse transcriptase. The PCR amplification process of complementary DNA can be monitored by fluorescent dyes. By detecting the total fluorescence of the product after PCR cycles, it can show whether the suspected specimen contains SARS-CoV-2 nucleic acid.

[1] WHO Timeline - COVID-19

<https://www.who.int/news-room/detail/27-04-2020-who-timeline---covid-19>

[2] COVID-19 Coronavirus Cases and Deaths https://www.worldometers.info/coronavirus/?utm_campaign=homeAdvegas17%22%20%5C%20%22countries%3Ca%20href=

[3] CDC Symptoms of COVID-19

<https://www.cdc.gov/coronavirus/2019-ncov/symptoms-testing/symptoms.html>


[4] Antibody responses to SARS-CoV-2 in patients of novel coronavirus disease 2019 <https://academic.oup.com/cid/advance-article/doi/10.1093/cid/ciaa344/5812996>

[5] CLSI M41-A, Viral Culture; Approved Guideline

<https://clsi.org/standards/products/microbiology/documents/m41/>

TEST PRINCIPLE

FORA COVID-19 Antigen Rapid Test is a lateral flow chromatographic immunoassay in a sandwich design with colloidal gold as an indicator. The FORA COVID-19 Antigen Rapid Test is designed to detect antigen from the SARS-CoV-2 in fresh nasopharyngeal swabs directly from patients who are suspected of COVID-19 by their healthcare provider. This test allows for the detection of SARS-CoV and SARS-CoV-2. The test detects, but does not differentiate, between the two viruses.



The test cassette consists of a test strip containing:

- 1) Conjugate pad : Anti SARS-CoV-2 N protein IgG CGC
- 2) Test line (T Line) : Anti-SARS-CoV-2 N protein IgG
- 3) Control line (C Line) : Control line antibody

* CGC: colloidal gold conjugation

Place a fresh nasopharyngeal swab specimen in an extraction tube with a correct volume of extraction buffer. During this time, the virus particles in the nasopharyngeal swab will be destroyed, and the internal viral nucleoprotein will be exposed. When a correct volume of extraction buffer with exposed viral nucleoprotein is dispensed into the sample well of the test cassette, the specimen migrates by capillary action along the test strip.

If present in the sample, SARS-CoV-2 nucleoprotein will bind Anti SARS-CoV-2 N protein IgG CGC and the complex will migrate to the T line.

The T line has Anti-SARS-CoV-2 N protein IgG fixed on the surface of the test strip. When the antibody and SARS-CoV-2 nucleocapsid protein association complex move to the T line, sandwich immune complexes are formed. Two antibodies sandwich SARS-CoV-2 nucleocapsid protein in the middle, and the colloidal gold particles are fixed on the T line. Aggregation occurs and produces a colored line, indicating a positive test result for SARS-CoV-2 nucleocapsid protein.

The C line has Control line antibody fixed on the surface of the test strip. The C Line is an internal control which should exhibit a colored line by immune complexes regardless of the color development on the T line. If no C Line is observed, the test result is invalid and the specimen must be re-collected and tested with a new cassette.

1) Accessories Included with the Test:

20-Test Kit/Box:

- Individually Foil Packaged Test Cassettes (20)
Each cassette includes:
 - 1) Conjugate pad: Anti SARS-CoV-2 N protein IgG CGC
 - 2) Test line (T Line): Anti-SARS-CoV-2 N protein IgG
 - 3) Control line (C Line): Control line antibody
- Extraction Buffer (6 ml per bottle) (2): Detergent, salts and non-reactive ingredients
- Extraction Tube (20)
- Nozzle Cap (20)
- Sterile Nasopharyngeal Swabs (20) (Optional)
- User Manual (1)
- Quick Reference Guide (1)

2) Accessories Required But Not Included with the Test:

- Timer or watch.
- Any necessary personal protective equipment.

WARNINGS AND PRECAUTIONS

1. For *in vitro* diagnostic use.
2. For prescription use only.
3. This test is intended for the detection of proteins from SARS-CoV-2, not for any other viruses or pathogens.
4. Sample collection and handling procedures require specific training and guidance.
5. When collecting a nasopharyngeal swab sample, use the sterile nasopharyngeal swab supplied in the kit.
6. Do not use the kit contents beyond the expiration date printed on the label.
7. Please check the package before use. Do not use if the package is damaged or the seal is broken.
8. Discard and do not use any damaged or dropped Test Cassette or material.
9. Wear suitable protective clothing, gloves, and eye/face protection when handling the contents of this kit.
10. Use of Nitrile, Latex (or equivalent) gloves is recommended when handling patient samples.
11. Always wear gloves to handle patient samples.
12. Avoid splashing and formation of droplets.
13. The user should never open the foil package of the Test Cassette exposing it to the ambient environment until the Test Cassette is ready for immediate use.
14. Start the assay procedure immediately after removing the cassette from the foil package.
15. Do not reuse the used Test Cassette, Extraction Tubes or Nasopharyngeal Swabs.
16. The Extraction Buffer contains a salt solution (saline). If the solution is in contact with the skin or eye, flush with copious amounts of water.
17. To obtain accurate results, the User Manual instructions must be followed. Incorrect sampling or procedure may result in inaccurate test results.
18. To obtain accurate results, do not use visually bloody or overly viscous samples.
19. Dispose of containers and unused contents in accordance with Local regulatory requirements.
20. Use appropriate disinfectants to thoroughly remove spills.
21. The possibility of infection cannot be totally ruled out. Therefore, all materials should be handled with care and treated like specimens. In the event of exposure, follow the instructions of local regulations.
22. Specimens, reagent kits and materials that may be contaminated during inspection are considered infectious waste, and must be discarded in accordance with local biological infectious regulations.
23. Do not interchange or mix different specimens.
24. Wash hands thoroughly after handling.
25. Use certified sterile rayon, foam, polyester or flocked flexible-shaft NP swabs to collect a nasopharyngeal sample.


STORAGE

1. Store at 2-30°C. Avoid direct sunlight.
2. Kit contents are stable until the expiration date printed on the label.
3. The test cassette must be kept in the sealed foil package. After unpacking, use the test cassette immediately.
4. Avoid freezing or heating test cassettes or kit contents.


SAMPLE COLLECTION AND HANDLING

Nasopharyngeal Swab Sample

Use only nasopharyngeal swab samples for testing.



1. Use the nasopharyngeal swab supplied in the kit. To collect a nasopharyngeal swab sample, carefully insert the swab into the nostril that presents the most secretion under visual inspection. Gently rotate while inserting the swab into the palate (not upwards) until resistance is encountered.



2. Gently rub and roll the swab over the surface of the posterior nasopharynx.
3. Slowly remove the swab while rotating it.

- For details of COVID-19 specimen collection, please refer to the interim guidelines issued by CDC: Interim Guidelines for Collecting, Handling, and Testing Clinical Specimens for COVID-19. <https://www.cdc.gov/coronavirus/2019-nCoV/lab/guidelines-clinical-specimens.html>


- It is recommended that the collection of specimens be done by physicians or healthcare providers.
- The training in specimen collection is highly recommended because of the importance of specimen quality. If the correct procedure in handling specimens is not followed, the value of the test results may be compromised or even negated.

- The test should be conducted as soon as possible after the sample is collected.
- For optimal test performance, use the swabs supplied in the kit.
- Please follow the local regulations for the collection, storage and transport of the specimens.


ASSAY PROCEDURE

- Bring the extraction buffer and individually foil packaged test cassette to room temperature approximately 30 minutes before performing the assay.
- All specimens, assay materials and procedures must be handled at room temperature.
- Check the expiration date printed on each kit contents before use. Do not use any accessory past the expiration date.
- Remove the cassette from the foil package before use and place it on a flat and dry surface.


1. Identify the Cassette for each sample with the individual's name and/or number.
2. Add 10 drops (about 500 µL) of extraction buffer to the extraction tube until the liquid level reaches the mark on the bottom of the extraction tube.





3. Immerse the patient nasopharyngeal swab sample into the extraction tube. Roll the swab three-five (3-5) times while pressing the head against the bottom and side of the extraction tube. Try to dissolve the sample in the extract as much as possible. Leave the swab in the extraction buffer for 1 minute.



4. When removing the swab, roll the swab head toward the inside of the extraction tube and squeeze the sides of the tube to extract the liquid from the swab. Dispose of the used swab in your biohazardous waste.
5. Insert the nozzle cap and press tightly.



6. Add 3 drops about 100 µL of the processed sample into the sample (S) well. Do not handle or move the cassette until the test is completed and ready for interpretation.
7. An interpretation is available within 15-20 minutes. Some positive results may appear sooner.




CAUTION:
Do not read the results after 20 minutes. It may provide false results.

INTERPRETATION OF RESULTS

Valid Assay:

Positive:

In addition to the presence of the colored C line, if the colored T line appears, the test result indicates the presence of SARS-CoV-2 virus in the nasopharyngeal swab sample. The result is COVID-19 positive or COVID-19 reactive. Within the specified observation time, a very weak colored line should be judged as a positive result. False positive results may occur due to cross-reacting antigens from previous infections, such as other coronaviruses, or from other causes. Samples with positive results should be confirmed with a molecular diagnostic test (e.g. RT-PCR) and clinical findings before a diagnostic determination is made.




Negative:

If only the colored C line appears, the test result indicates that SARS-CoV-2 virus is not detected at the time when the nasopharyngeal swab sample was collected. The result is COVID-19 negative or COVID-19 non-reactive. Negative results do not rule out SARS-CoV-2 infection, particularly for patients who have been in contact with known infected persons or in areas with high prevalence of active infection. Follow-up testing with a molecular diagnostic test (e.g. RT-PCR) is necessary to rule out infection in these individuals.



Invalid Assay:

There should always be a colored control line in the control region regardless of the test result. If the control line is not seen, repeat the assay with a new test cassette.



Negative / High CT value:

Within a specified observation time, a very weak color on the T line should be judged as negative or high CT value. Please refer to the RT-PCR result.



QUALITY CONTROL

FORA COVID-19 Antigen Rapid Test uses the Internal Control as the mechanism for quality control. A colored Control (C) line is an internal procedural control. It confirms sufficient sample volume, adequate membrane wicking and correct procedural technique.

External positive and negative controls are not supplied with this kit; however, external positive and negative controls should be tested in consistent with good laboratory practice to confirm the test procedure and to verify proper test performance.

LIMITATIONS OF THE PROCEDURE

- The contents in this kit are to be used for the qualitative detection of SARS-CoV-2 antigens from nasopharyngeal swab.
- Failure to follow the test procedure or incorrect interpretation of results may adversely affect test performance and result in invalid interpretation.
- This test detects both viable (live) and non-viable, SARS-CoV, and SARS-CoV-2. The test performance depends on the amount of virus (antigen) in the sample and may or may not correlate with viral culture results performed on the same sample.
- A negative test result may occur if the amount of virus (antigen) in a sample is below the limit of the assay or if the sample was not collected properly.
- Test results must be evaluated in conjunction with other clinical data available to the physician.
- The color of the test line has no correlation with clinical symptoms and severity. The interpretation of the test results must be evaluated together with epidemiology, clinical symptoms, and other diagnostic methods.
- Positive test results do not rule out co-infections with other viruses.
- Positive test results do not differentiate between SARS-CoV and SARS-CoV-2.
- Negative results cannot completely rule out the possibility of COVID-19 infection. The possible cause is that the amount of virus (antigen) in the sample is too low to be detected or the sample is not collected properly. Negative results must be determined with an FDA authorized molecular assay.
- Users should test samples as quickly as possible after sample collection.
- If the differentiation of specific SARS viruses and strains is needed, additional testing, in consultation with local public health departments, is required.

ANALYTICAL PERFORMANCE

Limit of Detection (LoD)

The Limit of Detection (LoD) of FORA COVID-19 Antigen Rapid Test was determined using limiting dilutions of live SARS-CoV-2, isolate TWN/CG-MH-CGU-01. The material was supplied frozen at a concentration of $10^{5.4}$ TCID₅₀ per mL. The study to determine the FORA COVID-19 Antigen Rapid Test LoD was designed to reflect the assay when using direct nasopharyngeal swab.

In this study, all the SARS-CoV-2 serial dilutions were made in the SARS-CoV-2 negative nasopharyngeal swab pool.

The LoD was determined in three steps:

- LoD Screening**
10-fold dilutions of the live SARS-CoV-2 were made as described above. These dilutions were tested in triplicate. The concentration demonstrating 3 of 3 positives was chosen for LoD range finding. Based on this testing, the concentration chosen for LoD Range Finding was $10^{2.4}$ TCID₅₀ per mL.
- LoD Range Finding**
Five (5) 2-fold dilutions of the $10^{2.4}$ TCID₅₀ per mL concentration were made as described above. These dilutions were tested in triplicate. The concentration demonstrating 3 of 3 positives was chosen for LoD confirmation. Based on this test the concentration chosen was 1.26×10^2 TCID₅₀ per mL.
- LoD Confirmation**
The concentration 1.26×10^2 TCID₅₀ per mL dilution was tested for a total of twenty (20) results. Twenty (20) of twenty (20) results were positive. Based on this test the concentration of LoD was confirmed as 1.26×10^2 TCID₅₀ per mL.

Cross-Reactivity

Cross-reactivity of the FORA COVID-19 Antigen Rapid Test was evaluated by testing various viruses (16) and bacteria (18). Each virus or bacteria was tested in triplicate in the absence or presence of 3.78×10^2 TCID₅₀/mL (3 LoD) of live SARS-CoV-2. The final concentration of each virus or bacteria was listed in the Table below. Testing was performed in triplicate.

Based on the data generated by this study, each virus or bacteria tested with FORA COVID-19 Antigen Rapid Test does not cross-react or interfere.

Clinical Performance

Clinical performance of FORA COVID-19 Antigen Rapid was determined by testing 57 positive and 80 negative specimens for SARS CoV-2 antigen (Ag) to have a sensitivity of 95.8% (95% CI: 86.0%-98.9%) and specificity of 98.6% (95% CI: 92.6%-100%).

FORA COVID-19 Antigen Rapid Test (TD-4531)	PCR Test Result			Subtotal
	Positive	Negative		
	Positive	53	1	
Negative	4 (2 subjects Ct>30)	79	83	
Subtotal	57	80	137	

Cross-Reactivity: FORA COVID-19 Antigen Rapid Test - Wet Testing

Virus/Bacteria	Concentration	Cross-Reactive Results	SARS-CoV-2 Concentration (3 LoD)	Interference Results
Human Coronavirus OC43	2.5 x 10 ⁵ pfu/mL	Negative	3.78 x 10 ² TCID ₅₀ /mL	Positive
Human Coronavirus 229E	2.5 x 10 ⁵ pfu/mL	Negative	3.78 x 10 ² TCID ₅₀ /mL	Positive
Influenza A + H1N1	2.5 x 10 ⁵ pfu/mL	Negative	3.78 x 10 ² TCID ₅₀ /mL	Positive
Influenza A + H3N2	2.5 x 10 ⁵ pfu/mL	Negative	3.78 x 10 ² TCID ₅₀ /mL	Positive
Influenza B + Victoria	2.5 x 10 ⁵ pfu/mL	Negative	3.78 x 10 ² TCID ₅₀ /mL	Positive
Influenza B + Yamagata	2.5 x 10 ⁵ pfu/mL	Negative	3.78 x 10 ² TCID ₅₀ /mL	Positive
Respiratory syncytial virus	2.5 x 10 ⁵ pfu/mL	Negative	3.78 x 10 ² TCID ₅₀ /mL	Positive
Rhinovirus	2.5 x 10 ⁵ pfu/mL	Negative	3.78 x 10 ² TCID ₅₀ /mL	Positive
Adenovirus type 1 (Adenoid 71)	2.5 x 10 ⁵ pfu/mL	Negative	3.78 x 10 ² TCID ₅₀ /mL	Positive
Adenovirus type 7	2.5 x 10 ⁵ pfu/mL	Negative	3.78 x 10 ² TCID ₅₀ /mL	Positive
Enterovirus 68	2.5 x 10 ⁵ pfu/mL	Negative	3.78 x 10 ² TCID ₅₀ /mL	Positive
Human parainfluenza type 1	2.5 x 10 ⁵ pfu/mL	Negative	3.78 x 10 ² TCID ₅₀ /mL	Positive
Human parainfluenza type 2	2.5 x 10 ⁵ pfu/mL	Negative	3.78 x 10 ² TCID ₅₀ /mL	Positive
Human parainfluenza type 3	2.5 x 10 ⁵ pfu/mL	Negative	3.78 x 10 ² TCID ₅₀ /mL	Positive
Human parainfluenza type 4	2.5 x 10 ⁵ pfu/mL	Negative	3.78 x 10 ² TCID ₅₀ /mL	Positive
Respiratory syncytial virus type B	2.5 x 10 ⁵ pfu/mL	Negative	3.78 x 10 ² TCID ₅₀ /mL	Positive
Bordetella pertussis	2 x 10 ⁶ CFU/mL	Negative	3.78 x 10 ² TCID ₅₀ /mL	Positive
Chlamydia pneumoniae	2 x 10 ⁶ CFU/mL	Negative	3.78 x 10 ² TCID ₅₀ /mL	Positive
Corynebacterium sp.	2 x 10 ⁶ CFU/mL	Negative	3.78 x 10 ² TCID ₅₀ /mL	Positive
Escherichia coli	2 x 10 ⁶ CFU/mL	Negative	3.78 x 10 ² TCID ₅₀ /mL	Positive
Hemophilus influenzae	2 x 10 ⁶ CFU/mL	Negative	3.78 x 10 ² TCID ₅₀ /mL	Positive
Lactobacillus sp.	2 x 10 ⁶ CFU/mL	Negative	3.78 x 10 ² TCID ₅₀ /mL	Positive
Moraxella catarrhalis	2 x 10 ⁶ CFU/mL	Negative	3.78 x 10 ² TCID ₅₀ /mL	Positive
Mycobacterium tuberculosis (avirulent)	2 x 10 ⁶ CFU/mL	Negative	3.78 x 10 ² TCID ₅₀ /mL	Positive
Neisseria meningitidis	2 x 10 ⁶ CFU/mL	Negative	3.78 x 10 ² TCID ₅₀ /mL	Positive
Neisseria sp.	2 x 10 ⁶ CFU/mL	Negative	3.78 x 10 ² TCID ₅₀ /mL	Positive
Pseudomonas aeruginosa	2 x 10 ⁶ CFU/mL	Negative	3.78 x 10 ² TCID ₅₀ /mL	Positive
Staphylococcus aureus (Protein A producer)	2 x 10 ⁶ CFU/mL	Negative	3.78 x 10 ² TCID ₅₀ /mL	Positive
Staphylococcus epidermidis	2 x 10 ⁶ CFU/mL	Negative	3.78 x 10 ² TCID ₅₀ /mL	Positive
Streptococcus pneumoniae	2 x 10 ⁶ CFU/mL	Negative	3.78 x 10 ² TCID ₅₀ /mL	Positive
Streptococcus pyogenes	2 x 10 ⁶ CFU/mL	Negative	3.78 x 10 ² TCID ₅₀ /mL	Positive
Streptococcus salivarius	2 x 10 ⁶ CFU/mL	Negative	3.78 x 10 ² TCID ₅₀ /mL	Positive
Pooled human nasal wash – representative of normal respiratory microbial flora	2 x 10 ⁶ CFU/mL	Negative	3.78 x 10 ² TCID ₅₀ /mL	Positive
Mycoplasma pneumoniae	2 x 10 ⁶ CFU/mL	Negative	3.78 x 10 ² TCID ₅₀ /mL	Positive

Interference Substances Studies

The study performed demonstrates that twenty (20) potentially interfering substances that may be found in the upper respiratory tract do not cross-react or interfere with the detection of SARS-CoV-2 in FORA COVID-19 Antigen Rapid Test. Each substance was tested in triplicate in the absence or presence of 3.78×10^2 TCID₅₀/mL (3 LoD) of live SARS-CoV-2.

Based on the data generated by this study, the substances tested in FORA COVID-19 Antigen Rapid Test do not cross-react or interfere.

Interfering Substance	Active Ingredient	Concentration	Cross-Reactive Results	SARS-CoV-2 Concentration (3 LoD)	Interference Results
Ephrine Nasal Spray "GPCPC"	Oxymetazoline	5% v/v	Negative	3.78 x 10 ² TCID ₅₀ /mL	Positive
Chloraseptic, Regular strength	Benzocaine / Menthol	1.5 mg/mL	Negative	3.78 x 10 ² TCID ₅₀ /mL	Positive
Tamiflu	Oseltamivir	2.5 mg/mL	Negative	3.78 x 10 ² TCID ₅₀ /mL	Positive
Physiomer Saline nasal spray	Saline	15% v/v	Negative	3.78 x 10 ² TCID ₅₀ /mL	Positive
Tobrex Eye Ointment	Tobramycin	51.4 µmol/L	Negative	3.78 x 10 ² TCID ₅₀ /mL	Positive
Sucrets	Dyclonine / Menthol	1.5 mg/mL	Negative	3.78 x 10 ² TCID ₅₀ /mL	Positive
NeilMed NasoGEL Spray	sodium hyaluronate / Saline	5% v/v	Negative	3.78 x 10 ² TCID ₅₀ /mL	Positive
Acetaminophen	Acetaminophen	1324 µmol/L	Negative	3.78 x 10 ² TCID ₅₀ /mL	Positive
Acetylsalicylic acid	Acetylsalicylic acid	3.62 mmol/L	Negative	3.78 x 10 ² TCID ₅₀ /mL	Positive
Ibuprofen	Ibuprofen	2.425 mmol/L	Negative	3.78 x 10 ² TCID ₅₀ /mL	Positive
Erythromycin	Erythromycin	81.6 µmol/L	Negative	3.78 x 10 ² TCID ₅₀ /mL	Positive
Fisherman's Friend	Menthol	1.5 mg/mL	Negative	3.78 x 10 ² TCID ₅₀ /mL	Positive
Plaquenil	Hydroxychloroquine sulphate	150 µmol/L	Negative	3.78 x 10 ² TCID ₅₀ /mL	Positive
SUPEROCIN	Ciprofloxacin	30.2 µmol/L	Negative	3.78 x 10 ² TCID ₅₀ /mL	Positive
Zeffix	Lamivudine	1 mg/mL	Negative	3.78 x 10 ² TCID ₅₀ /mL	Positive
Blood (human)	Blood (human)	2.5% v/v	Negative	3.78 x 10 ² TCID ₅₀ /mL	Positive
Ricola	Menthol	1.5 mg/mL	Negative	3.78 x 10 ² TCID ₅₀ /mL	Positive
Mupirocin	Mupirocin	10 mg/mL	Negative	3.78 x 10 ² TCID ₅₀ /mL	Positive
Flonase	Fluticasone	5% v/v	Negative	3.78 x 10 ² TCID ₅₀ /mL	Positive
Purified mucin protein	Mucin protein	2.5 mg/mL	Negative	3.78 x 10 ² TCID ₅₀ /mL	Positive

SYMBOL INFORMAION

SYMBOL	REFERENT	SYMBOL	REFERENT
	In vitro diagnostic medical device		Do not re-use
	Use-by date		Consult instructions for use
	Batch code		Manufacturer
	Temperature limit		Authorized representative in the European Community
	CE mark		Prescription use only
	Do not use if package is damaged		Model number
	Contains sufficient for <n> tests		

MedNet EC-REP GmbH
Borkstraße 10, 48163
Münster, Germany

TaiDoc Technology Corporation
B1-7F, No.127, Wugong 2nd Rd.,Wugu
Dist., 24888 New Taipei City, Taiwan

For Professional Use Only

