
# FORA®

# Diamond VOICE

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
 FORA DM20a / FORA DM20b  
For self-testing / Pour auto-mesure

Version 4.0 2021/08 311-4282100-008

# Dear FORA Diamond VOICE System Owner:

Thank you for purchasing the **FORA Diamond VOICE** Blood Glucose Monitoring System. This manual provides important information to help you use the system properly. Before using this product, please read the following contents thoroughly and carefully.

Regular monitoring of your blood glucose levels can help you and your doctor gain better control of your diabetes. Due to its compact size and easy operation, you can use the **FORA Diamond VOICE** Blood Glucose Monitoring System to easily monitor your blood glucose levels by yourself anywhere, any time.

If you have other questions regarding this product, please contact the place of purchase or call the Customer Care Line.

# **IMPORTANT SAFTY INSTRUCTIONS READ BEFORE USE**

1. Use this device **ONLY** for the intended use described in this manual.
2. Do **NOT** use accessories which are not specified by the manufacturer.
3. Do **NOT** use the device if it is not working properly or if it is damaged.
4. Do **NOT** use the equipment in places where aerosol sprays are being used or where oxygen is being administered.
5. Do **NOT** under any circumstances use the device on newborns or infants.
6. This device does **NOT** serve as a cure for any symptoms or diseases. The data measured is for reference only.
7. Before using this device to test blood glucose, read all instructions thoroughly and practice the test. Carry out all the quality control checks as directed.
8. Keep the device and testing equipment away from young children. Small items such as the battery cover, batteries, test strips, lancets and vial caps are choking hazards.
9. Use of this instrument in a dry environment, especially if synthetic materials are present (synthetic clothing, carpets etc.) may cause damaging static discharges that may cause erroneous results.
10. Do **NOT** use this instrument in close proximity to sources of strong electromagnetic radiation, as these may interfere with the accurate operation.

**KEEP THESE INSTRUCTIONS IN A SAFE PLACE**

# TABLE OF CONTENTS

<b>BEFORE YOU BEGIN</b>	<b>04</b>
Important Information	04
Intended Use	05
Test Principle	05
Meter Overview	06
Display Screen	07
Test Strip	08
Setting the Meter	09
<b>THE FOUR MEASURING MODES</b>	<b>13</b>
Control Solution Testing	14
Performing a Control Solution Test	14
<b>TESTING WITH BLOOD SAMPLE</b>	<b>17</b>
Preparing the Puncture Site	17
Alternative Site Testing	19
Performing a Blood Glucose Test	21
<b>METER MEMORY</b>	<b>24</b>
Reviewing Test Results	24
Reviewing Blood Glucose Day Average Results	25
Downloading Results onto a Computer	26
<b>MAINTENANCE</b>	<b>29</b>
Battery	29
Caring for Your Meter	30
Caring for Your Test Strips	31
Important Control Solution Information	31
<b>SYSTEM TROUBLESHOOTING</b>	<b>32</b>
Result Readings	32
Error Messages	33
Troubleshooting	35
<b>DETAILED INFORMATION</b>	<b>36</b>
<b>SYMBOL INFORMATION</b>	<b>37</b>
<b>SPECIFICATIONS</b>	<b>38</b>

# BEFORE YOU BEGIN

## Important Information

- Severe dehydration and excessive water loss may cause readings which are lower than actual values. If you believe you are suffering from severe dehydration, consult a healthcare professional immediately.
- If your blood glucose results are lower or higher than usual, and you do not have any symptoms of illness, first repeat the test. If you have symptoms or continue to get results which are higher or lower than usual, follow the treatment advice of your healthcare professional.
- Use only fresh whole blood samples to test your blood glucose. Using other substances will lead to incorrect results.
- If you are experiencing symptoms that are inconsistent with your blood glucose test results and you have followed all the instructions given in this owner's manual, contact your healthcare professional.
- We do not recommend using this product on severely hypotensive individuals or patients in shock. Readings which are lower than actual values may occur for individuals experiencing a hyperglycaemic-hyperosmolar state, with or without ketosis. Please consult the healthcare professional before use.
- The measurement unit used for indicating the concentration of blood or plasma glucose can either have a weight dimension (mg/dL) or a molarity (mmol/L). The approximate calculation rule for conversion of mg/dL in mmol/L is:

mg/dL	Divided by 18	= mmol/L
mmol/L	Times 18	= mg/dL

For example:

1)  $120 \text{ mg/dL} \div 18 = 6.6 \text{ mmol/L}$

2)  $7.2 \text{ mmol/L} \times 18 = 129 \text{ mg/dL}$  approximately.

## Intended Use


This system is intended for use outside the body (*in vitro* diagnostic use) by people with diabetes at home and by health care professionals in clinical settings as an aid to monitoring the effectiveness of diabetes control. It is intended to be used for the quantitative measurement of glucose (sugar) in samples of venous whole blood and fresh capillary whole blood (from the finger, palm, forearm or upper arm). It should not be used for the diagnosis of diabetes, or testing on newborns.

Professionals may test with capillary and venous blood sample; home use is limited to capillary whole blood testing. Use only heparin for anticoagulation of whole blood.

## Test Principle

Your system measures the amount of sugar (glucose) in whole blood. The glucose testing is based on the measurement of electrical current generated by the reaction of glucose with the reagent of the strip. The meter measures the current, calculates the blood glucose level, and displays the result. The strength of the current produced by the reaction depends on the amount of glucose in the blood sample.

## Meter Overview

**1 Test Slot**

Insert test strip here to turn the meter on for testing.

**2 Test Strip Ejector**

Eject the used strip by pushing up this button.

**3 Main Button**

Enter the meter memory and silence a reminder alarm.

**4 Display Screen**

**5 Data Port/USB (for DM20a)**

Download test results with a cable connection.

**6 Bluetooth Switch (for DM20b)**

Activating the Bluetooth function.

**7 Bluetooth Indicator (for DM20b)**


Download test results with a Bluetooth connection.

**8 SET Button**

Enter and confirm the meter settings.

## Display Screen

- |   |                       |    |                      |
|---|-----------------------|----|----------------------|
| 1 | Volume Symbol         | 9  | Time                 |
| 2 | Test Result           | 10 | Day Average          |
| 3 | Ketone Warning        | 11 | Alarm Symbol         |
| 4 | Warning Symbol        | 12 | Face/Low/High Symbol |
| 5 | Testing Mode Symbol   | 13 | Measurement Unit     |
| 6 | Memory Symbol         | 14 | Low Battery Symbol   |
| 7 | Control Solution Mode | 15 | Test Strip Symbol    |
| 8 | Date                  | 16 | Blood Drop Symbol    |

## Test Strip

### Absorbent Hole

Apply a drop of blood here. The blood will be automatically absorbed.

### Confirmation Window

This is where you confirm if enough blood has been applied to the absorbent hole in the strip.

### Test Strip Handle

Hold this part to insert the test strip into the slot.

### Contact Bars

Insert this end of the test strip into the meter. Push it in firmly until it will go no further.

### ATTENTION:

**The front side of test strip should face up when inserting test strip.**

Test results might be wrong if the contact bar is not fully inserted into the test slot.

### NOTE:

The **FORA Diamond VOICE** meter should only be used with **FORA Diamond VOICE** Test Strips. Using other test strips with this meter can produce inaccurate results.


## Setting the Meter

Before your first time using the meter and after every time changing the battery, you should check and update these settings. Make sure you complete the steps below and have your desired settings saved.

### Entering the Setting Mode

Start with the meter off (no test strip inserted). **Press SET.**

#### 1. Setting the date

With the year flashing, press Main button until the correct year appears. Press SET.


With the month flashing, press Main button until the correct month appears. Press SET.

With the day flashing, press Main button until the correct day appears. Press SET.

#### 2. Setting the time format

Press Main button to select the desired time format — 12h or 24h. Press SET.

### 3. Setting the time

With the hour flashing, press Main button until the correct hour appears. Press SET.



With the minute flashing, press Main button until the correct minute appears. Press SET.

### 4. Deleting the memory

With “dEL” and a flashing “**M**” on the display, press Main button and select “no” to keep the results in memory then press SET to skip.

To delete all the results, press Main button and “yes” and “**M**” are displayed on the meter, press SET to delete the memory.



### 5. Setting the speaking volume

There are seven (7) speaking volume options to choose from. Press Main button until the desired speaking volume appears. To confirm your selection, press SET.

Volume 0 indicates that the speaking function is turned off, and “” will not display during testing. Volume 1 to 7 indicates speaking volume from low to

high, and “🔊” will be displayed throughout the testing.



## 6. Choosing a language

Press the Main button to select L1/L2/L3. The default language for the meter is L2, which is French. To confirm your selection, press SET.

L1: English


L2: French (default setting)

L3: Arabic

## 7. Setting the reminder alarm

You may set up any or all of the reminder alarms (1-4). The meter displays “On” or “OFF” and “🔔”, press Main button to turn on or turn off to set the first reminder alarm.

Press Main button to select “On”, then press SET to set the hour. When the hour is flashing, press Main button to add an hour. Press SET to confirm and go to minutes, press Main button to add one minute. Hold Main button longer to add faster. Press SET to confirm and go to the next alarm setting.

If you do not want to set an alarm, press SET to skip this step.

If you want to turn off an alarm, find the alarm number by pressing SET in the setting mode, press Main button to change from “ON” to “OFF”.

At the time of your alarm, the meter will beep and automatically turn on. You can press Main button to silence the alarm and insert a test strip to begin testing. If you do not press Main button, the meter will beep for 2 minutes then switch off. If you do not want to test at this time, press Main button to switch off the meter.

### 8. Setting the Auto-Send (for DM20b)

Press Main button to select the auto-send On or OFF. Press SET.

**NOTE:** This function is referring to the Bluetooth data transmission. If “On” is selected, your result will be transmitted automatically right after the test.

**Congratulations! You have completed all settings!**

**NOTE:**


- These parameters can **ONLY be changed** in the setting mode.
- If the meter is idle for 3 minutes during the setting mode, it will switch off automatically.

# THE FOUR MEASURING MODES

The meter provides you with four modes for measuring, General, AC, PC and QC.

MODES	USE WHEN
General (displays as “Gen”)	any time of day without regard to time since last meal
AC	no food intake for at least 8 hours
PC	2 hours after a meal
QC	testing with the control solution

You can switch between each mode by:



1. Start with the meter switched off. Insert a test strip to turn on the meter. The screen will display a flashing “” and “Gen”.

## Voice Guide

Please apply blood onto the strip.

2. Press Main button to switch between General, AC, PC and QC mode.

## Control Solution Testing

Our Control Solution contains a known amount of glucose that reacts with test strips and is used to ensure your meter and test strips are working together correctly.

Test strips, control solutions, or sterile lancets may not be included in the kit (please check the contents on your product box). They can be purchased separately.

### Do a control solution test when:

- you suspect the meter or test strips are not working properly,
- your blood glucose test results are not consistent with how you feel, or if you think the results are not accurate,
- you have dropped or think you may have damaged the meter.

## Performing a Control Solution Test

You will need:


**1 Meter**

**2 Test Strip**





**3 Control Solution**

### 1. Insert the test strip to turn on the meter

Insert the test strip into the meter. Wait for the meter to display “” and “.



## 2. Press Main button to mark this test as a control solution test

With “**QC**” displayed, the meter will not store your test result in memory. If you press the Main button three times, the “**QC**” will disappear and this test is no longer a control solution test.


### WARNING:

When doing the control solution test, you have to mark it so that the test result will **NOT** mix with the blood glucose **TEST RESULTS** stored in the memory. Failure to do so will mix up the blood glucose test results with the control solution test results in memory.

## 3. Apply control solution

Shake the control solution vial thoroughly before use. Squeeze out a drop and wipe it off, then squeeze out another drop and place it on the tip of the vial cap. Hold the meter to move the absorbent hole of the test strip to touch the drop. Once the confirmation window fills completely, the meter will begin counting down. To avoid contaminating the control solution, do not directly apply control solution onto a strip.

#### 4. Read and compare the result

After counting down to 0, the control solution test result will appear on the display. Compare this result with the range printed on the test strip vial and it should fall within this range. If not, please read the instructions again and repeat the control solution test.

(112 mg/dL = 6.2 mmol/L; 100–165 mg/dL = 5.6–9.2 mmol/L)

#### Voice Guide

Your blood glucose is (number) milligram per deciliter.

#### Out-of-range results

If you continue to have test results fall outside the range printed on the test strip vial, the meter and strips may not be working properly. Do **NOT** test your blood. Contact the local customer service or place of purchase for help.

#### NOTE:

- The control solution range printed on the test strip vial is for control solution use only. It is not a recommended range for your blood glucose level.
- See the **Meter Maintenance** section for important information about your control solutions.

# TESTING WITH BLOOD SAMPLE

## **WARNING:**

To reduce the chance of infection:

- Never share a lancet or the lancing device.
- Always use a new, sterile lancet. Lancets are for single use only.
- Avoid getting hand lotion, oils, dirt, or debris in or on the lancets and the lancing device.

## **Preparing the Lancing Device for Blood Testing**

Please follow the instructions in the lancing device insert for collecting a blood sample.

## **Preparing the Puncture Site**

Stimulating blood perfusion by rubbing the puncture site before blood extraction has a significant influence on the glucose value obtained. Blood from a site that has not been rubbed exhibits a measurably different glucose concentration than blood from the finger. When the puncture site was rubbed prior to blood extraction, the difference was significantly reduced.

### **Please follow the suggestions below before obtaining a drop of blood:**

- Wash and dry your hands before starting.
- Select the puncture site either at fingertips or another body parts (please see section “Alternative Site Testing” (AST) on how to select the appropriate sites).
- Clean the puncture site using cotton moistened with 70% alcohol and **let it air dry.**
- Rub the puncture site for about 20 seconds before penetration.
- Use a clear cap (included in the kit) while setting up the lancing device.


### ● **Fingertip testing**

Press the lancing device's tip firmly against the lower side of your fingertip. Press the release button to prick your finger, then a click indicates that the puncture is complete.

### ● **Blood from sites other than the fingertip**

Replace the lancing device cap with the clear cap for AST. Pull the cocking control back until it clicks. When lancing the forearm, upper arm or hand, avoid lancing the areas with obvious veins because of excessive bleeding.


#### **NOTE:**

- Choose a different spot each time you test. Repeated punctures at the same spot may cause soreness and calluses.
- Please consult your health care professional before you begin AST.
- It is recommended that you discard the first drop of blood as it might contain tissue fluid, which may affect the test result.

## Alternative Site Testing

Alternative site testing (AST) is when individuals check their blood glucose levels using other areas of the body other than the fingertips. The FORA test strips allow AST to be performed on sites other than the fingertips.

### **Important:**

**There are limitations with AST (Alternative Site Testing).**

**Please consult your health care professional before you perform AST.**

### **When to use AST?**

Alternative site sample results may be different from fingertip sample results when glucose levels are changing rapidly (e.g., after a meal, after taking insulin, or during or after exercise).

We strongly recommend that you perform AST **ONLY** at the following times:

- In a pre-meal or fasting state (more than 2 hours since the last meal).
- Two hours or more after taking insulin.
- Two hours or more after exercise.

Do **NOT** rely on test results at an alternative sampling site, but use samples taken from the fingertip, if any of the following applies:

- You think your blood sugar is low.
- You are not aware of symptoms when you become hypoglycemic.
- The results do not agree with the way you feel.
- After a meal.
- After exercise.
- During illness.
- During times of stress.

To obtain a blood sample from the alternative sites, please rub the puncture site for approximately 20 seconds.

1. Replace the lancing device cap with the clear cap.
2. Pull the cocking control out until the orange bar appears on the release button window.

### **Important!**

- Do not use results from alternative site samples to calibrate continuous glucose monitoring systems (CGMS), or for insulin dose calculations.
- Choose a different spot each time you test. Repeated punctures at the same spot may cause soreness and calluses.
- Avoid lancing the areas with obvious veins to avoid excessive bleeding.
- It is recommended to discard the first drop of blood as it might contain tissue fluid, which may affect the test result.

# Performing a Blood Glucose Test

You will need:

**1 Meter**


**2 Test Strip**

**3 Lancing Device**

**4 Lancet**

## 1. Insert the test strip to turn on the meter


Wait for the meter to display “” and “”.

### Voice Guide

Thank you for using this product. Please relax during the measurement. Please apply blood onto the test strip.

## 2. Select the appropriate measuring mode by pressing Main button


For selecting the measurement mode, please refer to the “FOUR MEASURING MODES”.

## 3. Obtaining a blood sample

Use the pre-set lancing device to puncture the desired site. Wipe off the first appeared drop of blood with a clean cotton swab. The size of the drop should be at least as big as , which is **0.5** microliter ( $\mu\text{L}$ ) of volume. Gently squeeze the punctured area to obtain

another drop of blood. Be careful **NOT** to smear the blood sample.


#### 4. Apply the sample

Gently apply the drop of blood to the absorbent hole of the test strip at a titled angle. Confirmation window should be completely filled if enough blood sample has been applied. Do **NOT** remove your finger until you hear a beep sound.

#### Voice Guide

Please apply blood onto the strip.

#### NOTE:

- Do not press the punctured site against the test strip or try to smear the blood.
- If you do not apply a blood sample to the test strip within 3 minutes, the meter will automatically turn off. You must remove and reinsert the test strip to start a new test.
- The confirmation window should be filled with blood before the meter begins to count down. **NEVER** try to add more blood to the test strip after the drop of blood has moved away. **Discard the used test strip and retest with a new one.**
- If you have trouble filling the confirmation window, please contact your health care professional or the local customer service for assistance.


## 5. Read your result

The result of your blood glucose test will appear after the meter counts down to 0. The blood glucose result will be stored in the memory automatically.

(100 mg/dL = 5.6 mmol/L)

### Voice Guide

The blood glucose is (number) milligrams per deciliter.

## 6. Eject the used test strip

Eject the test strip by pushing the eject button on the side. Use a sharp bin to dispose of used test strips. The meter will switch itself off automatically.

**Always follow the instructions in the lancing device insert when removing the lancet.**


### **WARNING:**

The used lancet and test strip may be biohazardous . Please discard them carefully according to your local regulations.

# METER MEMORY

The meter stores the 450 most recent blood glucose test results along with respective dates and times in its memory. To enter the meter memory, **start with the meter switched off.**


## Reviewing Test Results

### 1. Press and release Main button

“M” will appear on the display. Press Main button again, and the first reading you see is the last blood glucose result along with date, time and the measurement mode.

(300 mg/dL = 16.7 mmol/L)


### 2. Press the Main button to recall the test results stored in the meter each time you press. After the last test result, press Main button again and the meter will be turned off.

(200 mg/dL = 11.1 mmol/L)

### Voice Guide

Your blood glucose is (number) milligrams per deciliter.

## Reviewing Blood Glucose Day Average Results

### 1. Press and release Main button.

When “**M**” appears on the display, keep pressing Main button for 3 seconds until the flashing “**DAY AVG**” appears. Release Main button and then your 7-day average result measured in general mode will appear on the display. (100 mg/dL = 5.6 mmol/L)

**2. Press Main button to review** 14-, 21-, 28-, 60- and 90- day average results stored in each measuring mode in the order of Gen, AC, then PC.

**3. Exit the meter memory.** Keep pressing the Main button and the meter will turn off after displaying the last test result.

### NOTE:

- Any time you wish to exit the memory, keep pressing Main button for 5 seconds or leave it without any action for 3 minutes. The meter will switch off automatically.
- Control solution results are **NOT** included in the day average.
- If using the meter for the first time, “—” displays when you recall the test results or review the average result. It indicates that there is no test result in the memory.

## Downloading Results onto a Computer


**FORA Diamond VOICE** has 2 types of transmission methods; your meter uses either USB or Bluetooth to transmit the data, please check your meter box for the transmission method of your meter.

### Data Transmission Via Cable (For DM20a)

You can use the meter with an USB cable and the Health Care Software System to view test results on your personal computer. To learn more about the Health Care Software System or to obtain an USB cable separately, please contact local customer services or the place of purchase for assistance.

#### 1. Obtaining the required cable and installing the software

To download the Health Care Software System, please visit the Foracare Suisse AG, website: [www.foracare.ch](http://www.foracare.ch)

#### 2. Connecting to a personal computer

Connect the cable to a USB port on your computer. With the meter switched off, connect the other end of the USB cable to the meter data port. “USB” will appear on the meter display, indicating that the meter is in communication mode.

#### 3. Data transmission

To transmit data, follow the instructions provided with the software. Results will be transmitted with date and time. Remove the cable and the meter will automatically switch off.

## Data Transmission Via Bluetooth (For DM20b)

You can use your device with an iOS (5.0.1 or higher), Android system (4.3 API Level 18 or higher) or PC (Windows 8 or higher) to download data from your meter via Bluetooth. Follow the steps below to transmit data from your meter. Please contact local customer service or place of purchase for assistance.

Please note that you must complete the pairing between meter and Bluetooth receiver before transmitting data.

1. Turn on the Bluetooth switch on the right side of your blood glucose meter and your mobile device.
2. Your mobile device will now begin searching for the Bluetooth signal to pair with.
3. Once the mobile device has found the blood glucose meter; the name of the meter will appear on the pairing list. Please select and add to the list.
4. On the mobile phone, the blood glucose meter will be shown as a paired device and on the blood glucose meter the Bluetooth Indicator LED will lit solid, indicating that the devices have been successfully paired.

### Bluetooth indicator on the blood glucose meter:

BLUETOOTH INDICATOR	STATUS
Flash Fast	The meter is searching the device of Bluetooth signal.
Flash Slowly	The meter is pairing with the device of Bluetooth.
Lit Solid	The meter is transmitting the data now the connection is established.

## **Mobile phone compatibility issues**

The Bluetooth functionality is implemented in different ways by the various mobile phone manufacturers.

Unfortunately, in some mobile phone models, even with Bluetooth functionality, they may be compatible only with certain types of devices. If a problem occurs in the connection between your mobile phone and the blood glucose meter, or if you are uncertain regarding your mobile phone's Bluetooth capabilities, please consult your mobile phone manual or blood glucose meter manual or alternatively contact your local customer service for assistance.

### **WARNING:**

- While the meter is in transmission mode, it will be unable to perform a blood glucose test.
- Make sure your device has turned on Bluetooth before transmitting the data and the meter is within the receiving range.

# MAINTENANCE

## Battery

Your meter comes with two 1.5V AAA size alkaline batteries.

### Low Battery Signal

The meter will display one of the messages below to alert you when the meter power is getting low.

The “” appears with **E-b, Error and low:**

The power is not enough to do a test. Please change the batteries immediately.

### Replacing the Battery

To replace the batteries, make sure that the meter is turned off.

1. Press the edge of the battery cover and lift it up to remove.
2. Remove the old batteries and replace with two 1.5V AAA size alkaline batteries.
3. Close the battery cover. If the batteries are inserted correctly, you will hear a “beep” afterwards.

#### NOTE:

- Replacing the batteries does not affect the test results stored in the memory.

- As with all small batteries, these batteries should be kept away from children. If swallowed, promptly seek medical assistance.
- Batteries might leak chemicals if unused for a long time. Remove the batteries if you are not going to use the device for an extended period (i.e., 3 months or more).
- Properly dispose of the batteries according to your local environmental regulations.

## Caring for Your Meter

To avoid the meter and test strips attracting dirt, dust or other contaminants, please wash and dry your hands thoroughly before use.

### Cleaning

- To clean the meter exterior, wipe it with a cloth moistened with tap water or a mild cleaning agent, then dry the device with a soft dry cloth. Do **NOT** rinse with water.
- Do **NOT** use organic solvents to clean the meter

### Meter Storage

- Storage conditions: -20 °C to 60 °C (-4 °F to 140 °F), below 95% relative humidity.
- Always store or transport the meter in its original storage case.
- Avoid dropping and heavy impact.
- Avoid direct sunlight and high humidity.

### Meter Disposal

- The used meter should be treated as contaminated and may carry a risk of infection during measurement. The batteries in this used meter should be removed and the meter should be disposed in accordance with local regulations.
- The meter falls outside the scope of the European Directive 2012/19/EU-Directive on waste electrical and electronic equipment (WEEE).

## Caring for Your Test Strips

- Storage conditions: 2°C to 32°C (35.6°F to 89.6°F), below 85% relative humidity. Do **NOT** freeze.
- Store your test strips in their original vial only. Do not transfer to another container.
- Store test strip packages in a cool dry place. Keep away from direct sunlight and heat.
- After removing a test strip from the vial, immediately close the vial cap tightly.
- Touch the test strip with clean and dry hands.
- Use each test strip immediately after removing it from the vial.
- Write the opening date on the vial label when you first opened it. Discard remaining test strips after 6 months.
- Do not use test strips beyond the expiry date. This may cause inaccurate results.
- Do not bend, cut, or alter a test strip in any way.
- Keep the strip vial away from children since the cap and the test strip may be a choking hazard. If swallowed, promptly see a doctor for help.

For further information, please refer to the test strip package insert.

## Important Control Solution Information

- Use only our control solutions with your meter.
- Do not use the control solution beyond the expiry date or 3 months after first opening. Write the opening date on the control solution vial and discard the remaining solution after 3 months.
- It is recommended that the control solution test be done at room temperature 20°C to 25°C (68°F to 77°F). Make sure your control solution, meter, and test strips are at this specified temperature range before testing.
- Shake the vial before use, discard the first drop of control solution, and wipe off the dispenser tip to ensure a pure sample and an accurate result.
- Store the control solution tightly closed at temperatures between 2°C to 30°C (36°F to 86°F). Do **NOT** freeze.

# SYSTEM TROUBLESHOOTING

If you follow the recommended action but the problem persists, or error messages other than the ones below appear, please call your local customer service. Do not attempt to repair yourself and never try to disassemble the meter under any circumstances.

## Result Readings

MESSAGE	WHAT IT MEANS		
<b>Lo</b>	< 20 mg/dL (1.1 mmol/L)		
<b>LOW</b> ☹️	20–69 mg/dL (1.1–3.8 mmol/L)		
😊	<b>AC</b> 🍏	<b>PC</b> 🍷	<b>Gen</b>
	70–129 mg/dL (3.9–7.2 mmol/L)	70–179 mg/dL (3.9–9.9 mmol/L)	70–119 mg/dL (3.9–6.6 mmol/L)
<b>HIGH</b> ☹️	<b>AC</b> 🍏	<b>PC</b> 🍷	<b>Gen</b>
	130–239 mg/dL (7.2–13.3 mmol/L)	180–239 mg/dL (9.9–13.3 mmol/L)	120–239 mg/dL (6.7–13.3 mmol/L)
<b>KETONE?</b> ☹️	≥ 240 mg/dL (13.3 mmol/L)		
<b>Hi</b>	> 600 mg/dL (33.3 mmol/L)		

## Error Messages

Message	What does meter say?	Cause	What To Do
	Only Error Beep	There is insufficient power in the batteries to perform the test.	Replace the batteries immediately and reset the date and time on the meter setting.
	Only Error Beep	The strip has been used.	Repeat the test with a new strip.
	Only Error Beep	The blood sample has been applied to the strip before the meter was ready for the measurement.	Repeat the test with a new strip. Apply the blood sample when the meter indicates that it is ready for the measurement by displaying a “  ” and a “  ” which starts blinking.
	Only Error Beep	Ambient temperature is below the operating temperature.	The operating temperature range is from 10 °C to 40 °C (50 °F to 104 °F). Repeat the test after the device and test strip have reached the above-mentioned operating temperature.
	Only Error Beep	Ambient temperature is above the operating temperature.	
	Only Error Beep	You may have removed the strip after applying blood, or there is insufficient blood volume.	Review the instructions and repeat the test with a new test strip. Please make sure you use the correct technique and dose with enough blood.

Message	What does meter say?	Cause	What To Do
	Only Error Beep	Meter calibration data error.	Review the instructions and repeat the test with a new strip and using the correct technique. If the problem persists, contact your local customer service for assistance.
	Only Error Beep	The component of the device might be broken.	
	None	Communicate to voice chip error.	
	Only Error Beep	Meter malfunction due to damaged components.	
	Only Error Beep	Incorrect data on the meter or the code strip.	

## Troubleshooting

1. If the meter does not display a message after inserting a test strip:

POSSIBLE CAUSE	WHAT TO DO
Battery exhausted.	Replace the battery.
Test strip inserted upside down or incompletely.	Insert the test strip with contact bars end first and facing up.
Defective meter or test strips.	Please contact customer services.

2. If the test does not start after applying the sample:

POSSIBLE CAUSE	WHAT TO DO
Insufficient blood sample.	Repeat the test using a new test strip with larger volume of blood sample.
Defective test strip.	Repeat the test with a new test strip.
Sample applied after automatic switch-off (2 minutes after last user action).	Repeat the test with a new test strip. Apply sample only when flashing “  ” appears on the display.
Defective meter.	Please contact customer services.

3. If the control solution testing result is out of range:

POSSIBLE CAUSE	WHAT TO DO
Error in performing the test.	Read instructions thoroughly and repeat the test again.
Control solution vial was poorly shaken.	Shake the control solution vigorously and repeat the test again.
Expired or contaminated control solution.	Check the expiry date of the control solution.
Control solution that is too warm or too cold.	Control solution, meter, and test strips should be at room temperature 20°C to 25°C (68°F to 77°F) before testing.
Defective test strip.	Repeat the test with a new test strip.
Meter malfunction.	Please contact customer services.

## DETAILED INFORMATION


The meter provides you with plasma equivalent results.

Time of day	Normal plasma glucose range for people <b>without</b> diabetes (mg/dL)
Fasting and before meal	< 100 mg/dL (5.6 mmol/L)
2 hours after meals	< 140 mg/dL (7.8 mmol/L)

Source: American Diabetes Association (2012). Clinical Practice Recommendations. Diabetes Care, 35 (Supplement 1): S1–S100.

**Please consult your doctor to determine a target range that works best for you.**

# SYMBOL INFORMATION

SYMBOL	REFERENT
	<i>In vitro</i> diagnostic medical device
	Do not reuse
	Consult instructions for use
	Keep away from sunlight
	Keep dry
	Temperature limitation
	Use by/ Expiry date
	Humidity limitation
	Use within 6 months after first opening
	Batch code
	Authorised representative in the European Union

SYMBOL	REFERENT
	Caution, consult accompanying documents
	Sterilized using irradiation
	Do not use if package is damaged
	Model No.
	Disposal of waste equipment
	Battery
	3 Volts DC
	Dispose of the packaging properly after use
	CE Mark
	Serial number
	Manufacturer
	Medical device

# SPECIFICATIONS

**Model No.:** FORA DM20a / FORA DM20b

**Dimension & Weight:** 104.4 x 58 x 23 mm, 71.2 g

**Power Source:** Two 1.5V AAA alkaline batteries

**Display:** LCD

**Memory:** 450 measurement results with respective date and time

**External Output:** USB / Bluetooth

Auto electrode insertion detection

Auto sample loading detection

Auto reaction time count-down

Auto switch-off after 3 minutes without action

Temperature warning

**Operating Condition:**

10 °C to 40 °C (50 °F to 104 °F), below 85% R.H. (non-condensing)

**Storage/Transportation Conditions:**

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**ForaCare Suisse AG**
**MedNet EC-REP GmbH**