

# FORA® IR21 / IR20

## IR21

Ear/Forehead Thermometer /  
Ohr-/Stirnthermometer /  
Thermomètre Auriculaire/Frontal /  
Termometro Auricolare/Frontale /  
Termómetro de frente y oído /  
Termómetro de Ouvido/Testa

## IR20

Ear Thermometer /  
Ohrthermometer /  
Thermomètre Auriculaire /  
Termometro Auricolare /  
Termómetro de oído /  
Termómetro de Ouvido

**Owner's Manual / Benutzerhandbuch /  
Manuel de l'utilisateur / Manuale del Proprietario /  
Manual de instrucciones / Manual de instruções**

# IMPORTANT SAFETY INSTRUCTIONS

English

## READ THIS BEFORE USING

The following basic safety precautions should always be taken.

1. Close supervision is necessary when the thermometer is used by, on, or near children, handicapped persons or invalids.
2. Use the thermometer only for the intended use described in this manual.
3. Do not use the thermometer if it is not working properly, or if it has suffered any damage.
4. Do not use accessories which are not supplied or recommended by the manufacturer.

**KEEP THESE INSTRUCTIONS IN A SAFE PLACE**

## Cautions and Warnings

- ▲ As with any thermometer, proper technique is crucial to getting accurate temperature readings. **Please read this manual thoroughly and carefully before using.**
- ▲ Always store the thermometer in a cool and dry place: temperatures between  $-20^{\circ}\text{C}$  to  $60^{\circ}\text{C}$  ( $-4^{\circ}\text{F}$  to  $140^{\circ}\text{F}$ ) relative humidity less than 95%. **Avoid direct sunlight.**
- ▲ **Avoid dropping** the thermometer from a height or strongly hitting it with a hard object.
- ▲ **Do not** touch the probe lens.
- ▲ **Do not** disassemble the thermometer.
- ▲ Basic safety precautions should always be observed, especially when the thermometer is used on or near children and disabled persons.
- ▲ This thermometer is not intended to be a substitution for consultation with your physician.
- ▲ The skin/surface scan temperature serves as a reference only. It cannot be a judgment on fever.
- ▲ Temperature of left and right ear may differ. Always measure using same ear.

- ▲ Only use the probe cover for infrared FORA IR20 ear thermometer.
- ▲ For proper hygiene, do not share probe cover. Damaged probe cover may result in error display.
- ▲ If you experience any serious incident that occurred in relation to the use of this product, please report it to the manufacturer and the competent authority of medical devices in your country.
  
- ▲ A serious incident means any incident that directly or indirectly led, might have led, or might lead to any of the following:
  - (a) the death of a patient, user, or other people,
  - (b) the temporary or permanent serious deterioration of a patient's, user's or other person's state of health,
  - (c) a serious public health threat.

## Restrictions of Use

This thermometer is clinically proven to produce accurate temperature measurements. However, please be advised if you have the following situations:

- ▲ The accuracy cannot be ensured for a person who has deformity in the ear such that the thermometer probe

can not be properly inserted into the ear canal.

- ▲ The accuracy cannot be ensured when blood or drainage is found in the ear canal.
- ▲ Take temperature from the other ear if ear drops or medications have been placed in an ear.
- ▲ For a person who wears ear plug or hearing aid, remove the device and wait for 15 minutes before taking temperature.

### **NOTE**

*Never try to clean inside the ears. You may accidentally damage the eardrum or its surrounding tissues. Remove excess earwax only when you can reach it with a clean cloth. Consult a physician if you suspect the presence of excess earwax.*

## Introduction

Thank you for choosing FORA IR21 ear/forehead thermometer. This innovative medical device relies on advanced infrared (IR) technology to measure temperature instantly and accurately through two convenient sites: the ear canal and the surface of the skin over the temple. FORA IR21 ear/forehead thermometer is an elegantly designed infrared thermometer meant for your family.

Thank you for choosing FORA IR20 ear thermometer. This innovative medical device relies on advanced infrared (IR) technology to measure temperature instantly. FORA IR20 ear thermometer is an elegantly designed infrared thermometer meant for your family.

## Intended use

FORA IR21 ear/forehead thermometer is intended for the intermittent measurement and monitoring of human body temperature from ear canal or skin/surface scan. The device is intended for use by people of all ages in the home.

FORA IR20 ear thermometer is intended for the intermittent measurement and monitoring of human body temperature from ear canal. The device is intended for use by people of all ages in the home.


## How does it work

The thermometer measures the infrared heat generated by the eardrum and its surrounding tissue, or by the surface of the skin over the temporal artery. The thermometer then converts it into a temperature value shown on LCD.


### NOTE

*The thermometer does not emit any infrared signal.*

## Thermometer Parts

## LCD Screen


- ① Ear temperature indicator
- ② Skin/Surface temperature indicator
- ③ Temperature scanning in progress
- ④ Temperature display
- ⑤ Memory mode
- ⑥ Record numbers
- ⑦ Low battery indicator
- ⑧ Temperature unit
- ⑨ Communication symbol
- ⑩ Non-probe cover warning (FORA IR20)

## Replacing the Battery

The thermometer comes with two 1.5V AAA alkaline batteries. Replace it when “” appears. Please follow the steps to replace new batteries.


### Step 1

Remove the battery cover.


### Step 2

Place the new battery in the battery compartment and press it in until the battery is firmly secured.

## Step 3

Reattach the battery cover.

### NOTE

- ▲ Although the thermometer works when “” appears, we still recommend that you change the battery to obtain an accurate result.
- ▲ Remove the battery if stored for a long period of time.
- ▲ The battery should be kept out of reach of children. If they are swallowed, promptly see a doctor for help.

## About Normal Body Temperature & Fever

Body temperature can vary from one individual/person to next. It also varies by location on the body and different time of the day. Below shows the statistical normal ranges from different sites. Please keep in mind that temperatures measured from different sites, even at the same time, should not be directly compared.

Fever indicates that the body temperature is higher than normal. This symptom may be potentially caused by infection, overdressing or immunization. Some people may not experience fever even when they are ill. These include, but are not limited to, infants younger than 3 months old, persons with compromised immune systems, persons taking antibiotics, steroids or antipyretics (aspirin, ibuprofen, acetaminophen), or persons with certain chronic illnesses. **Please consult your physician when you feel ill even if you do not have fever.**

**Table 1<sup>1</sup> Body Site Normal Temperature Range**

Oral	0.6°C (1°F) or more above or below 37°C (98.6°F)
Rectal/ear	0.3°C to 0.6°C (0.5°F to 1°F) higher than oral temperature
Axillary (underarm)	0.3°C to 0.6°C (0.5°F to 1°F) lower than oral temperature

\*1. Body Temperature at WebMD;

<http://firstaid.webmd.com/body-temperature> retrieved at 2010 Jan 7.

## Using the Device

### As an Ear Thermometer (IR21)

## Step 1


Remove the probe cover first.

## Step 2

Turn on the thermometer.

Press and release the On / Memory button. When ready, the thermometer displays the last measurement.

## Step 3


Gently fit the probe into the ear canal.

## Step 4

Press and release the Scan button.

Do not remove the thermometer until it beeps.

## Step 5

Read the result.

"" is shown together with a temperature value.

### NOTE

- ▲ The thermometer does not require any costly lens filters.
- ▲ Turn off the thermometer by pressing On / Memory button twice. It will automatically turn off if left idle for 3 minutes.
- ▲ If you need to take another reading, wait "" flashing before taking another measurement.

## Hints on Taking Ear Temperature

As with other of thermometer, you may observe slight variations in consecutive measurements. ***It is recommended that you take 3 temperature readings and use the highest one for the following situations:***

- ▲ Infants younger than 3 months old.
- ▲ Children younger than 3 years old and who have a compromised immune system and the presence / absence of fever is critical.
- ▲ When you are learning to use the thermometer.

**Do not take a reading while eating and /or talking. Wait 30 minutes after any of the following situations before taking a measurement:**

- ▲ When you have your ear covered.
- ▲ After exercising, swimming, or taking a bath. When expose to extreme temperature.

To take accurate readings, the ear must be free from excess earwax build-up.

## Replacing Probe Cover (IR20)

### Step 1

Check the probe cover is clean and new.

### Step 2

Attach new probe cover into the probe.

***Now, you can start to test your ear temperature.***

### NOTE

- ▲ *If a probe cover is not attached on the probe, LCD will show " ~~X~~ " until a new probe cover is firmly attached to the probe.*
- ▲ *The used probe cover should be removed after each measurement by a clean and new probe cover to ensure accuracy reading.*
- ▲ *For hygiene use, attach a new probe cover each time and do not touch its tip.*

## Taking Ear Temperature (IR20)


### Step 1

Make sure probe cover is firmly attached to the probe.

### Step 2

Press and release the On / Memory button to turn on the thermometer. When ready, the thermometer displays the last measurement.

### Step 3

Stretch your ear canal by pulling your ear backwards and upwards and carefully insert the thermometer into the ear.

### Step 4

Press and release the Scan button. Do not remove the thermometer until it beeps.

### Step 5

Read the result. "ⓘ" is shown together with a temperature value.

### Step 6

Discard used probe cover into trash by pressing probe cover ejection button.

### NOTE

Turn off the thermometer by pressing On / Memory button twice. It will automatically turn off if left idle for 3 minutes.

## As a Skin/Surface Scan Thermometer (IR21)

### Step 1

**Make sure the probe cover is capped.**

### Step 2

**Turn on the thermometer.**




Press and release the On / Memory button. The thermometer displays the last measurement.

### Step 3

**Put the probe on the temple.**

Press and hold the Scan button. Make sure the probe is flat on the temple, not at an angle.

### Step 4

**Read the result.**

Release the button and read the result. “☺” is shown together with a temperature value and beeps.

### NOTE


- ▲ As the skin/surface scan temperature is likely to be affected by sweat, oil and the surroundings, the reading shall be taken as a **reference only**.
- ▲ If the probe is placed at an angle on the skin/surface scan, the reading will be affected by surrounding temperature.
- ▲ You can perform a skin/surface scan from 3 cm away.
- ▲ Turn off the thermometer by pressing On / Memory button twice. Or leave for 3 minutes after a temperature measurement, the thermometer will automatically turn off.
- ▲ If you need to take another reading, for the flashing “▶▶▶” flashing before taking another measurement.

## Recalling the Memory

Your thermometer stores 10 recent readings in the memory.


### Step 1

**Be sure the thermometer is OFF before recalling this memory.**


### Step 2

**Press the On / Memory button to turn on the thermometer.**

### Step 3

**Press the On / Memory button for 3 seconds to enter memory mode.**

Each time you press the On / Memory button, a result will be displayed in the order of dates (latest result shown first), together with “M” and number ( from 1 to 10 ).

When the memory is full, the oldest result is deleted as the new one added. When the last record is displayed in the LCD, press On / Memory button again to return the first record.

### Step 4

**Exit the memory.**

Press the Scan button and LCD will show latest results with flashing “▶▶”. Then press On/Memory button again and LCD will show “OFF” to exit memory.

#### NOTE

*When using the thermometer for the first time, the first memory result will display 0.0°C. It indicates that there are no test results in memory.*

## Viewing Results on a Personal Computer

Results in memory can be transmitted to the personal computer by either cable or wireless connection

for FORA IR21c / IR20c and FORA IR21b / IR20b, respectively. Accessories needed to activate this function are:

- ▲ Health Care Software System: a software downloaded from ForaCare's website ([www.foracare.ch](http://www.foracare.ch))
- ▲ Interface Cable: an optional accessory for FORA IR21c / IR20c.
- ▲ Bluetooth Adapter: an optional accessory for FORA IR21b / IR20b.

Please contact your local customer service for above accessories.

## **Transmitting data via Cable (FORA IR21c / IR20c)**

### **Step 1 Install Software**

Install Health Care System Software on your computer by following the instructions provided on ForaCare's website.

### **Step 2 Connect to Personal Computer**

Connect the interface cable to a serial port on the back of your computer. With the thermometer turn off, connect

the interface cable to the data port located at the bottom of the thermometer. Then " PCL " will appear in the display, indicating that the thermometer is ready to transmit data.

## **Step 3 Transmit Data**

Follow the instructions provided in the software to transmit data. Results transmitted will include date and time. Remove the cable and the thermometer will automatically turn off.

## **Transmitting data via Bluetooth Adapter (FORA IR21b / IR20b)**

### **Step 1**

Install Health Care System Software on your computer by following the instructions provided on ForaCare's website.

### **Step 2**

Install Bluetooth Adapter Software.

### **Step 3**

Reset FORA IR21b / IR20b thermometer.

### **Step 4**


Pair FORA IR21b / IR20b thermometer with your computer. Select FORA IR21b / IR20b thermometer device code: "FORA-Device" and key in the Bluetooth passkey (PIN): 111111 for pairing.

## Step 5

Connect Bluetooth Serial Port Service:

Right click on icon "FORA-Device" and select "connect" and "Bluetooth Serial Port Service" for connection. The results will be transmitted automatically to your computer once you finish taking the temperature. In the meanwhile, the communication symbol on the thermometer will flash.

If you are using the Bluetooth Adapter provided by ForaCare Suisse AG, please follow the installation guide "Bluetooth Easy Card" within the package.

Bluetooth Indicator


Bluetooth Indicator	Description
Flash fast	The meter is searching the device of Bluetooth signal.
Flash slowly	The meter is pairing with the device of Bluetooth.
Lit solid	The meter is transmitting the data now the connection is completed.

## NOTE

- ▲ While the thermometer is connected to the PC or the server, it is unable to perform a test.
- ▲ We recommend that you use bluetooth adaptor provided by ForaCare Suisse AG.

## Care & Cleaning

- ▲ Keep the probe clean, as earwax and grease build up may affect the measurement.
- ▲ The probe tip is the most delicate part of the thermometer. It has to be clean and intact to ensure accurate readings. Very gently wipe the surface with a cotton swab.
- ▲ The body of the thermometer is not water-resistant. Never put the thermometer under a running tap or submerge it into water. Use a soft and dry cloth to clean it. Do not use abrasive cleaners.
- ▲ Store the thermometer in a cool and dry location. Free from dust and away from direct sunlight.











## Fever Indicator

When the measurement is equal or over 38°C (100.4°F), red backlight is shown together with the result. The red backlight will last for 3 seconds and turn off.

## Troubleshooting

The table below shows problems you may encounter. All error messages below would be shown together with orange backlight. Please follow “what to do” to resolve problems. If the problem still exists, please call your local dealer for help.

Message	What it means	What to do
	Room temperature is below 10°C.	Put the thermometer under operating temperature range of 10°C to 40°C (50°F to 104°F).
	Room temperature is above 40°C.	Put the thermometer under operating temperature range of 10°C to 40°C (50°F to 104°F).
 	You don't use probe cover while measuring ear temperature.	Please place probe cover into probe again.

   	Problem with the thermometer.	Review the instructions and re-start the measurement procedure. If the above steps do not work. Please contact the dealer.
	Battery is low and “  ” appears on LCD.	Please replace batteries as soon as possible.
	Appears when the batteries can't provide enough power for a test.	Please replace new batteries.
 	Temperature measurement falls outside the displayed temperature range: (ear temperature range from 32°C to 43°C. skin/surface scan temperature range from 22°C to 44°C.)	Please follow this manual to take a reading again.

# Specifications

Dimensions	162.2mm (L) x 32.4mm (W) x 38.0mm (H) (FORA IR21) 160.5mm (L) x 32.5mm (W) x 45.9mm (H) (FORA IR20)
Weight	110.0g (include 2 x1.5V AAA batteries) (FORA IR21) 105.5g (include 2 x1.5V AAA batteries) (FORA IR20)
Battery	2 x 1.5V AAA alkaline batteries
External output	RS-232 / USB / Bluetooth
Displayed temperature range	32°C to 43°C (89.6°F to 109.4°F)
Display resolution	0.1°C / 0.1°F
Accuracy	Meet the accuracy requirement specified in ASTM E1965-98 ▲ ±0.2°C (±0.4°F) for the range of 36.0°C to 39.0°C (96.8°F to 102.2°F) ▲ ±0.3°C (±0.5°F) from 34.0°C to 35.9°C (93.2°F to 96.6°F) and from 39.1°C to 42.2°C (102.4°F to 108.0°F)
Temperature unit	°C / °F
Operating temperature range	10°C to 40°C (50°F to 104°F)
Operating humidity	95% RH or less

Storage temperature range	-20°C to 60°C (-4°F to 140°F)
Storage humidity	95% RH or less
Memory capacity	10 measurements

The specifications may be changed without prior notice.

This device has been tested to meet the electrical and safety requirements of: ASTM E1965-98, EN 12470-5, EN 60601-1-2:2007/AC:2010, EN 60601-1-4:1996, EN 60601-1-6:2010, EN 301 489-17, EN 300 328.

# Symbol Information

Symbol	Referent
	Consult instructions for use
	Manufacturer
	Serial number
	Caution, consult accompanying documents
	Type BF Equipment
	CE Mark
	Disposal of waste equipment

	Storage/Transportation temperature limitation
	Storage/ Transportation humidity limitation
	Authorised representative in the European Union
	Medical Device

## Meter Disposal

The used meter should be treated as contaminated and may carry a risk of infection during measurement. The batteries in this used meter should be removed and the meter should be disposed in accordance with local regulations.

The meter falls outside the scope of the European Directive 2002/96/ EC-Directive on waste electrical and electronic equipment (WEEE).

## WARRANTY TERMS AND CONDITIONS

With respect to disposable products, ForaCare Suisse warrants to the original purchaser that, at time of delivery, each standard product manufactured by ForaCare Suisse shall be free of defects in material and workmanship and, when used for the purposes and indications described on the labeling, is fit for the purposes and indications described on the labeling. All warranties for a product shall expire as of the product expiration date, or if none, after two (2) year from the original date of purchase, as long as it has not been modified, altered, or misused. ForaCare Suisse warranty hereunder shall not apply if:

(i) a product is not used in accordance with its instructions or if it is used for a purpose not indicated on the labeling; (ii) any repairs, alterations or other work has been performed by the buyer or others on such item, other than work performed with ForaCare Suisse's authorisation and according to its approved procedures; or (iii) the alleged defect is a result of abuse, misuse, improper maintenance, accident or the negligence of any party other than ForaCare Suisse. The warranty set forth herein is conditioned upon proper storage, installation, use and maintenance in accordance with applicable written recommendations of ForaCare Suisse. The warranty furnished hereunder does not extend to damage items purchased hereunder resulting in whole or in part from the use of components, accessories, parts or supplies not furnished by ForaCare Suisse.

**Warning:** Medical electrical equipment needs special precautions regarding EMC and needs to be installed according to the EMC information provided. Careful consideration of this information is essential when stacking or collocating equipment and when routing cables and accessories.

**Warning:** RF mobile communications equipment can affect medical electrical equipment.

Recommended separation distance between portable and mobile RF communications equipment and FORA IR20/IR21			
FORA IR20/IR21 is intended for use in an electromagnetic environment (for home healthcare and professional healthcare) in which radiated RF disturbances are controlled. The customer or the user of FORA IR20/IR21 can help prevent electromagnetic interference by maintaining a minimum distance between the portable and mobile RF communications equipment (transmitters) and FORA IR20/IR21 as recommended below, depending on the maximum output power of the communications equipment.			
Rated maximum output power of transmitter W	Separation distance according to frequency of transmitter m		
	150 kHz to 80 MHz $d = 1,2\sqrt{P}$	80 MHz to 800 MHz $d = 1,2\sqrt{P}$	800 MHz to 2,7 GHz $d = 2,3\sqrt{P}$
0,01	N/A	0,12	0,23
0,1	N/A	0,38	0,73
1	N/A	1,2	2,3
10	N/A	3,8	7,3
100	N/A	12	23
For transmitters rated at a maximum output power not listed above, the recommended separation distance $d$ in meters (m) can be estimated using the equation applicable to the frequency of the transmitter, where $p$ is the maximum output power rating of the transmitter in watts (W) depending on the transmitter manufacturer.			
<b>NOTE 1</b> At 80 MHz and 800 MHz, the separation distance for the higher frequency range applies.			
<b>NOTE 2</b> These guidelines may not apply to all situations. Electromagnetic propagation is affected by absorption and reflection from structures, objects and people.			

**Manufacturer's declaration-electromagnetic emissions**

FORA IR20/IR21 is intended for use in the electromagnetic environment (for home healthcare and professional healthcare) specified below. The customer or the user of FORA IR20/IR21 should assure that it is used in such an environment.

<b>Emission test</b>	<b>Compliance</b>	<b>Electromagnetic environment-guidance</b> (for home healthcare and professional healthcare)
RF emissions CISPR 11	Group 1	FORA IR20/IR21 uses RF energy only for internal use. Therefore, its RF emissions are very low and are not likely to cause any interference from nearby electronic equipment.
RF emissions CISPR 11	Class B	FORA IR20/IR21 is suitable for use in all establishments, including domestic establishments and those directly connected to the public low-voltage power supply network that supplies buildings used for domestic purposes.
Harmonic emissions IEC 61000-3-2	Not applicable	
Voltage fluctuations/flicker emissions IEC 61000-3-3	Not applicable	

**Manufacturer's declaration-electromagnetic immunity**

FORA IR20/IR21 is intended for use in the electromagnetic environment (for home healthcare and professional healthcare) specified below.  
The customer or the user of FORA IR20/IR21 should assure that it is used in the environment specified below.

<b>Immunity test</b>	<b>IEC 60601 test level</b>	<b>Compliance level</b>	<b>Electromagnetic environment-guidance</b> (for home healthcare and professional healthcare environment)
Electrostatic discharge (ESD) IEC 61000-4-2	Contact: $\pm 8$ kV Air $\pm 2$ kV, $\pm 4$ kV, $\pm 8$ kV, $\pm 15$ kV	Contact: $\pm 8$ kV Air $\pm 2$ kV, $\pm 4$ kV, $\pm 8$ kV, $\pm 15$ kV	Floors should be wood, concrete or ceramic tile. If floors are covered with synthetic material, the relative humidity should be at least 30 %.
Electrical fast transient/burst IEC 61000-4-4	$\pm 2$ kV for power supply lines $\pm 1$ kV for input/output lines	Not applicable Not applicable	Mains power quality should be that of a typical home and professional healthcare environment.
Surge IEC 61000-4-5	$\pm 0.5$ kV, $\pm 1$ kV line(s) to line(s) $\pm 0.5$ kV, $\pm 1$ kV, $\pm 2$ kV line(s) to earth	Not applicable Not applicable	Mains power quality should be that of a typical home and professional healthcare environment.
Voltage Dips, short interruptions and voltage variations on power supply input lines IEC 61000-4-11	Voltage dips: 0 % $U_T$ ; 0,5 cycle 0 % $U_T$ ; 1 cycle 70 % $U_T$ ; 25/30 cycles  Voltage interruptions: 0 % $U_T$ ; 250/300 cycle	Voltage dips: Not applicable Not applicable Not applicable  Voltage interruptions: Not applicable	Mains power quality should be that of a typical home healthcare and professional healthcare environment. If the user of FORA IR20/IR21 requires continued operation during power mains interruptions, it is recommended that FORA IR20/IR21 be powered from an uninterruptible power supply or a battery.
Power frequency (50, 60 Hz) magnetic field IEC 61000-4-8	30 A/m 50 Hz or 60 Hz	30 A/m 50 Hz and 60 Hz	FORA IR20/IR21 power frequency magnetic fields should be at levels characteristic of a typical location in a typical home healthcare and professional healthcare environment.

**NOTE**  $U_T$  is the a.c. mains voltage prior to application of the test level.

**Manufacturer's declaration-electromagnetic immunity**

FORA IR20/IR21 is intended for use in the electromagnetic environment (for home healthcare and professional healthcare) specified below. The customer or the user of FORA IR20/IR21 should assure that it is used in the environment specified below.

Immunity test	IEC 60601 test level	Compliance level	Electromagnetic environment-guidance (for home healthcare and professional healthcare environment)
Conducted RF IEC 61000-4-6	3 Vrms: 0,15 MHz – 80 MHz 6 Vrms: in ISM and amateur radio bands between 0,15 MHz and 80 MHz  80 % AM at 1 kHz	Not applicable  Not applicable	<p><b>Portable and mobile RF communications equipment must not be used close to any parts of FORA IR20/IR21 including cables, other than the recommended separation distance calculated from the equation applicable to the frequency of the transmitter.</b></p> <p><b>Recommended separation distance:</b>  <math>d = 1,2 \sqrt{P}</math>  <math>d = 1,2 \sqrt{P}</math> 80MHz to 800 MHz  <math>d = 2,3 \sqrt{P}</math> 800MHz to 2,7 GHz</p> <p>Where <math>P</math> is the maximum output power rating of the transmitter in watts (<math>W</math>) according to the transmitter manufacturer and <math>d</math> is the recommended separation distance in metres (<math>m</math>).</p> <p>Field strengths from fixed RF transmitters, as determined by an electromagnetic site survey, <sup>a)</sup> should be less than the compliance level in each frequency range.<sup>b)</sup></p> <p>Interference may occur in the vicinity of equipment marked with the following symbol:</p> 
Radiated RF IEC 61000-4-3	10 V/m 80 MHz – 2,7 GHz 80 % AM at 1 kHz	10 V/m 80 MHz – 2,7 GHz 80 % AM at 1 kHz	

**NOTE 1** At 80 MHz and 800 MHz, the higher frequency range applies.

**NOTE 2** These guidelines may not apply to all situations. Electromagnetic propagation is affected by absorption and reflection from structures, objects and people.

**a)** Field strengths from fixed transmitters, such as base stations for radio (cellular/cordless) telephones and land mobile radios, amateur radio, AM and FM radio broadcast and TV broadcast cannot be predicted theoretically with accuracy. To assess the electromagnetic environment due to fixed RF transmitters, an electromagnetic site survey should be considered. If the measured field strength in the location in which FORA IR20/IR21 is used exceeds the applicable RF compliance level above, FORA IR20/IR21 should be observed to verify normal operation. If abnormal performance is observed, additional measures may be necessary, such as re-orienting or relocating FORA IR20/IR21.

**b)** Over the frequency range 150 kHz to 80 MHz, field strengths should be less than 3 V/m.

**Manufacturer's declaration-electromagnetic immunity**  
**Test specifications for ENCLOSURE PORT IMMUNITY to RF wireless communications equipment**

FORA IR20/IR21 is intended for use in the electromagnetic environment (for home healthcare and professional healthcare) specified below.  
 The customer or the user of FORA IR20/IR21 should assure that it is used in such an environment.

<b>Test frequency (MHz)</b>	<b>Band <sup>a)</sup> (MHz)</b>	<b>Service <sup>a)</sup></b>	<b>Modulation <sup>b)</sup></b>	<b>Maximum power (W)</b>	<b>Distance (m)</b>	<b>IMMUNITY TEST LEVEL (V/m)</b>	<b>Compliance LEVEL (V/m) (for home and professional healthcare)</b>
385	380 – 390	TETRA 400	Pulse modulation <sup>b)</sup> 18 Hz	1,8	0,3	27	27
450	430 – 47	GMRS 460, FRS 460	FM <sup>c)</sup> ±5 kHz deviation 1 kHz sine	2	0,3	28	28
710	704 – 787	LTE Band 13, 17	Pulse modulation <sup>b)</sup> 217 Hz	0,2	0,3	9	9
745							
780							
810	800 – 960	GSM 800/900, TETRA 800, iDEN 820, CDMA 850, LTE Band 5	Pulse modulation <sup>b)</sup> 18 Hz	2	0,3	28	28
870							
930							
1 720	1 700 – 1 990	GSM 1800; CDMA 1900; GSM 1900; DECT; LTE Band 1, 3, 4, 25; UMTS	Pulse modulation <sup>b)</sup> 217 Hz	2	0,3	28	28
1 845							
1 970							
2 450	2 400 – 2 570	Bluetooth, WLAN, 802.11 b/g/n, RFID 2450, LTE Band 7	Pulse modulation <sup>b)</sup> 217 Hz	2	0,3	28	28

5 240	5 100 – 5 800	WLAN 802.11 a/n	Pulse modulation <sup>b)</sup> 217 Hz	0,2	0,3	9	9
5 500							
5 785							

**NOTE** To achieve the IMMUNITY TEST LEVEL, the distance between the transmitting antenna and the ME EQUIPMENT or ME SYSTEM may be reduced to 1 m. The 1 m test distance is permitted by IEC 61000-4-3.

**a)** For some services, only the uplink frequencies are included.

**b)** The carrier shall be modulated using a 50 % duty cycle square wave signal.

**c)** As an alternative to FM modulation, 50 % pulse modulation at 18 Hz may be used because while it does not represent actual modulation, it would be worst case.

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**MedNet EC-REP GmbH**