FORA[®] COMFORT plus G30a

Blood Glucose Monitoring System
Blutzucker-Überwachungssystem
Système de contrôle de la glycémie
Sistema per il monitoraggio della glicemia

Owner's Manual
Bedienungsanleitung
Manuel de l'utilisateur
Manuale dell'utente

Safety Information

Read the following *Safety Information* thoroughly before using the device.

- Use this device **ONLY** for the intended use described in this manual.
- Do **NOT** use accessories which are not specified by the manufacturer.
- Do **NOT** use the device if it is not working properly or damaged.
- This device does **NOT** serve as a cure for any symptoms or diseases. The data measured is for reference only. Always consult your doctor to have the results interpreted.
- The blood glucose test strip can **NOT** be used for the testing of newborns.
- Before using this device to test blood glucose, read all instructions thoroughly and practice the test. Carry out all the quality control checks as directed.
- Keep the device and testing supplies away from young children. Small items such as the battery cover, batteries, test strips, lancets and vial caps are choking hazards.
- The use of this instrument in a dry environment, especially if synthetic materials are present (synthetic clothing, carpets etc.) may cause damaging static discharges that may cause erroneous results.
- Do **NOT** use this instrument in close proximity to sources of strong electromagnetic radiation, as these may interfere with the correct operation.

- Proper maintenance as well as timely calibration of the device together with the control solution is essential in ensuring the longevity of your device. If you are concerned about the accuracy of the measurement, please contact the place of purchase or customer service representative for assistance.

KEEP THESE INSTRUCTIONS IN A SAFE PLACE

Important Information

- Severe dehydration and excessive water loss may cause readings which are lower than actual values. If you believe you are suffering from severe dehydration, consult a healthcare professional immediately.
- If your blood glucose results are lower or higher than usual, and you do not have symptoms of illness, first repeat the test. If you have symptoms or continue to get results higher or lower than usual, follow the treatment advice of your healthcare professional.
- Use only fresh whole blood sample to test your blood glucose. Using other substances will lead to incorrect results.
- If you are experiencing symptoms that are inconsistent with your blood glucose test results and you have followed all instructions described in this owner's manual, call your healthcare professional.

- We do not recommend using this product on severely hypotensive individuals or patients in shock. Readings which are lower than actual values may occur for individuals experiencing a hyperglycaemic-hyperosmolar state, with or without ketosis. Please consult the healthcare professional before use.
- The measurement unit used for indicating the concentration of blood glucose can have mg/dL or mmol/L. The approximate calculation rule for conversion of mg/dL in mmol/L is:

| | | |
|--------|---------------|----------|
| mg/dL | Divided by 18 | = mmol/L |
| mmol/L | Times 18 | = mg/dL |

For example:

1) $120 \text{ mg/dL} \div 18 = 6.6 \text{ mmol/L}$

2) $7.2 \text{ mmol/L} \times 18 = 129 \text{ mg/dL}$ approximately

Introduction

Intended Use


This system is intended for use outside the body (in vitro diagnostic use) by people with diabetes at home and by health-care professionals in clinical settings as an aid to monitor the effectiveness of diabetes control. It is intended to be used for the quantitative measurement of glucose (sugar) in fresh capillary (from the finger, palm, forearm, and upper arm) blood samples. It should not be used for the diagnosis of diabetes, or testing on newborns.

Use only heparin for anticoagulation of whole blood.

Test Principle


Your system measures the amount of sugar (glucose) in whole blood. The glucose testing is based on the measurement of electrical current generated by the reaction of glucose with the reagent of the strip. The meter measures the current, calculates the blood glucose level, and displays the result. The strength of the current produced by the reaction depends on the amount of glucose in the blood sample.

Product Overview

1. Test Slot
2. Display Screen
3. **M** Button
4. **C** Button
5. Data Port
6. Test Strip Ejector
7. **SET** Button

Screen Display

- | | |
|-------------------------|----------------------|
| 1. Low Battery Symbol | 8. Measuring Mode |
| 2. Test Strip Symbol | 9. Date |
| 3. Blood Drop Symbol | 10. Time |
| 4. Test Result | 11. Ketone Warning |
| 5. Face/Low/High Symbol | 12. Error Message |
| 6. Memory Symbol | 13. Day Average |
| 7. Alarm Symbol | 14. Measurement Unit |


Getting Started

Initial Setup

Before using your meter for the first time or if you change the meter battery, you should check and update these settings. Make sure you complete the steps below and have your desired settings saved.

► Entering the Setting Mode.

Start with the meter off (no test strip inserted). Press **SET** to turn on the meter.


Step 1: Setting the date.

With the year flashing, press **M** until the correct year appears. Press **SET**.


With the month flashing, press **M** until the correct month appears. Press **SET**.

With the day flashing, press **M** until the correct day appears. Press **SET**.

Step 2: Setting the time format.


Press **M** and release to select the desired time format --- 12h or 24h. Press **SET**.

Step 3: Setting the time.

With the hour flashing, press **M** until the correct hour appears. Press **SET**.

With the minute flashing, press **M** until the correct minute appears. Press **SET**.


Step 4: Setting the Universal Tone

You may set up Beep On, Universal Tone On and Beep Off by press **M** to change it.

When Universal Tone is turned on, the meter guides you through the blood glucose test using beep tones; it also outputs the result as a series of beeps.

When Beep is turned off, the alarm function will remain effective.


Step 5: Deleting the memory.

With "dEL" and "M" symbol on the display, press **M** and select "no" to keep the results in memory then press **SET** to skip.

To delete all the results, press **M** and select "YES".


Then press **SET** button to delete all memory. "OK" is displayed in the meter, which indicates that all data stored is deleted.


Step 6: Setting the reminder alarm.

You may set up any or all of the reminder alarms (1-4). The meter displays "On" or "OFF" and "🕒", press **M** to turn on or turn off to set the first reminder alarm.

Press **M** to select "On", then press **SET** to set the hour. When the hour is flashing, press **M** to add an hour. Press **SET** to confirm and go to minutes, press **M** to add one minute. Hold **M** longer to add faster. Press **SET** to confirm and go to the next alarm setting.

If you do not want to set an alarm, press **SET** to skip this step.

If you want to turn off an alarm, find the alarm number by pressing **SET** in the setting mode, press **M** to change from "ON" to "OFF".

At the time of your alarm, the meter will beep and automatically turn on. You can press **M** to silence the alarm and insert a test strip to begin testing. If you do not press **M**, the meter will beep for 2 minutes then switch off. If you do not want to test at this time, press **M** to switch off the meter.

Press **SET** to turn off the meter. The meter repeats all your

settings before shutdown.

Congratulations! You have completed all settings!

Note:

- These parameters can **ONLY** be changed in the setting mode.
- If the meter is idle for 3 minutes during the setting mode, it will switch off automatically.

Before Testing

Control Solution Testing

Our Control Solution contains a known amount of glucose that reacts with test strips and is used to ensure your device and test strips are working together correctly.

Test strips, control solutions, or sterile lancets may not be included in the kit (please check the contents on your product box). They can be purchased separately. Please make sure you have those items needed for a blood glucose test beforehand.

Do a control solution test when:


- ✓ you first receive the device.
- ✓ you begin using a new vial of test strips.
- ✓ you suspect the device or test strips are not working properly.
- ✓ your blood glucose test results are not consistent with

- how you feel, or if you think the results are not accurate.
- ✓ you have dropped or think you may have damaged the device.

To perform the control solution test, do the following:


Step 1: Insert the test strip to turn on the meter.

Insert the test strip into the meter. Wait for the meter to display the test strip and blood drop symbols.

Step 2: Press C to mark this test as a control solution test.

If you press **C** again, the "QC" will disappear and this test is no longer a control solution test.

Step 3: Apply control solution.

Shake the control solution vial thoroughly before use. Squeeze out a drop and wipe it off, then squeeze out another drop and place it on the tip of the vial cap. Hold the meter to move the absorbent hole of the test strip to touch the drop. Once the confirmation window fills completely, the meter will begin counting down. To avoid contaminating the control solution do not directly apply control solution onto a strip.

Step 4: Read and compare the result.

After counting down to 0, the control solution test result will appear on the display. Compare this result with the range printed on the test strip vial and it should fall within this range. If not, please read the instructions again and repeat the control solution test.

Note:



- Do **NOT** test your blood.
- The control solution range printed on the test strip vial or individual foil pack is for control solution use only. It is not a recommended range for your blood glucose level.
- Refer to the **Maintenance** section for important information about your control solutions.

Out-of-range results:

If you continue to get results that fall outside the range printed on the test strip vial, it means that the meter and strips might not be working properly.

Testing Your Blood Glucose

Test Strip Appearance

ATTENTION:

The front side of test strip with FORA Logo should face up when inserting test strip. Test results might be wrong if the contact bar is not **fully** inserted into the test slot.

Inserting a Test Strip

Insert the test strip into its slot.

Important!

The front side of test strip should face up when inserting the test strip. Test results might be wrong if the contact bar is not fully inserted into the test slot.


Important!

To reduce the chance of infection:

- Never share a lancet or a lancing device.
- Always use a new, sterile lancet. Lancets are for single use only. Avoid getting hand lotion, oils, dirt, or debris in or on the lancets and the lancing device.

Preparing the Lancing Device


1. Remove the cap.

2. Insert a new lancet firmly into the white lancet holder cup.

3. Remove the protective disk on the lancet.

Hold the lancet firmly in place and twist off the protective disk.

4. Replace the cap until it snaps or clicks into place.

5. Rotate the dial to set the desired lancing depth.

6. Pull the cocking control out until the orange bar appears on the release button window.

Obtaining the Blood Sample


Please follow the suggestions below before obtaining a drop of blood:

- Wash and dry your hands before starting.
- Select the puncture site either on your fingertips or other body parts.

- Rub the puncture site for about 20 seconds before penetration.

❖ Blood from the fingertip

1. Press the lancing device tip firmly against the lower side of your fingertip.
2. Press the release button to prick your finger. A click indicates that the puncture is complete.

❖ Blood from sites other than the fingertip

Alternative site testing (AST) is when individuals check their blood glucose levels using other areas of the body other than the fingertips. The FORA test strips allow AST to be

performed on sites other than the fingertips. Please consult your health care professional before you begin AST.

We strongly recommend that you perform AST **ONLY** at the following times:


- During a pre-meal or fasting state (more than 2 hours since the last meal).
- Two hours or more after taking insulin.
- Two hours or more after exercise.

Do **NOT** use AST if:

- You think your blood glucose is low.
- You may not notice if you are hypoglycemic.
- Your AST results are inconsistent with the way you feel.
- You are testing for hyperglycemia.
- Your routine glucose results often fluctuate.

To obtain a blood sample from the alternative sites, please rub the puncture site for approximately 20 seconds.

1. Replace the lancing device cap with the clear cap.

2. Pull the cocking control out until the orange bar appears on the release button window.

Important!

- Choose a different spot each time you test. Repeated punctures at the same spot may cause soreness and calluses.
- Avoid lancing the areas with obvious veins to avoid excessive bleeding.
- It is recommended to discard the first drop of blood as it might contain tissue fluid, which may affect the test result.

Make Blood Glucose Test

1. Insert the test strip into the test slot of the device. Wait for the device to display the blood drop “”.

2. Select the appropriate measuring mode by pressing **M**.

- General Tests - any time of day without regard to time since the last meal.
- AC (**AC** ) - no food intake for at least 8 hours.
- PC (**PC** ) - 2 hours after a meal.
- QC (**QC**) - testing with control solution.


3. Obtain a blood sample.

Use the pre-set lancing device to puncture your desired site. After penetration, discard the first drop of blood with a clean tissue or cotton. Gently squeeze the punctured area to obtain another drop of blood. Be careful **NOT** to smear the blood sample. The volume of blood sample must be at least 0.5 microliter (μL) for blood glucose testing.


4. Apply the blood sample.

Move your finger to meet the absorbent hole of the test strip and the drop will be automatically be drawn onto the test strip. Remove your finger until the confirmation window is filled. The meter begins to count down. Do not remove your finger until you hear a beep sound.

5. Read your result.

The results of your blood glucose test will appear after the meter counts down to 0. The results will be stored automatically in the meter memory.

Announcement of Test Result by Universal Tone

The blood glucose results will be broken down into individual digits and each digit represents the corresponding number of beeps.

The result is announced three times in succession and each time is preceded by two quick beeps. So you will hear: 2 quick beeps – results – 2 quick beeps – results – 2 quick beeps – result.

For mg/dL meters, the hundreds are always announced, even when the result is below 100.

Examples:

80 mg/dL is announced as 1 long beep (0) – 1 single pause – 8 single beeps (8) – 1 single pause – 1 long beep (0)

182 mg/dL is announced as 1 single beep (1) – 1 single pause – 8 single beeps (8) – 1 single pause – 2 single beeps (2)

For mmol/L meters, the tens are always announced, even when the result is below 10. The decimal point is represented

by 1 quick beep.

Examples:

6.0 mmol/L is announced as 1 long beep (0) – 1 single pause – 6 single beeps (6) – 1 single pause – 1 quick beep (.) – 1 single pause – 1 long beep (0)

Note:

Information or warnings in the form of symbols displayed together with the results are not announced acoustically.

Eject the Used Test Strip and Remove the Lancet

To eject the test strip, point the strip at a disposal container for sharp items. The meter will switch itself off automatically after the test strip is ejected.

Important!

The used lancet and test strip may be biohazards. Please consult your health-care provider for proper disposal which complies with your local regulations.


Reviewing Test Results

Your device stores the 450 most recent test results along with respective dates and times in its memory. To enter the device memory, start with the device switched off.


To review all test results, do the following:

Step 1: Press M.

"**M**" will appear on the display. Press **M** again, and the first reading you see is the last blood glucose result along with date, time and the measuring mode.


Step 2: Press M to recall the test results stored in the meter each time you press.

Step 3: Exit the meter memory.


After the last test results, press **M** again and the meter will be turned off.



To review the day-average test results, do the following:

Step 1: Press M.

When "**M**" appears on the display, keep pressing **M** for 3 seconds until the flashing "**DAY AVG**" appears. Release **M** and then your 7-day average result measured in general mode will appear on the display.


Step 2: Press M to review 14-, 21-, 28-, 60- and 90- day average results stored in each measuring mode in the order of Gen, AC, then PC.


Step 3: Exit the meter memory.

Keep pressing **M** and the meter will turn off after displaying the last test result.

Note:

- Any time you wish to exit the memory, keep pressing **M** for 5 seconds or leave it without any action for 3 minutes. The meter will switch off automatically.
- Control solution results are NOT included in the day average.
- If using the meter for the first time, "----" displays when you recall the test results or review the average result. It indicates that there is no test result in the memory.

Announcement of Memory Result by Universal Tone

Only the most recent result that was saved can be announced acoustically. If you press **M** to turn the meter on, you will first hear the Long-BEEP which stands for power-on and then the most recent result.

Only the average for the last 7 days is announced acoustically. If the 7-day average cannot be calculated, three horizontal bars are displayed. This is signaled acoustically with 3 long beeps representing 3 zeroes.

Transferring Data

You can use the meter with an interface cable and the Health Care Software System to view your test results on your personal computer.

To learn more about the Health Care Software System or to obtain an interface cable separately, please contact local customer services or the place of purchase for assistance.


Data Transmission Via RS-232 Cable

Step 1: Obtaining the required cable and installing the software.

To download the Health Care Software System, please visit the Foracare Suisse AG, website at www.foracare.ch.

Step 2: Connecting to a personal computer.

Connect the cable to a serial port on your computer. With the meter switched off, connect the other end of the interface cable to the meter data port. "PC" will appear on the meter display, indicating that the meter is in communication mode.

Step 3: Data transmission.

To transmit data, follow the instructions provided with the software. Results will be transmitted with date and time. Remove the cable and the meter will automatically switch off.

WARNING!

While the meter is connecting to the PC, it will be unable to perform a blood glucose test.

Maintenance

Battery


Your meter comes with one 3V CR2032 lithium battery.

► Low Battery Signal

The "🔋" symbol appears along with display messages: The meter is functional and the result remains accurate, but it is time to change the battery.

The "🔋" symbol appears with E-b, Error and low: The power is not enough to do a test. Please change the battery immediately, automatically be switched off.

Low Battery Warning:

When the batteries will soon be exhausted, 2 quick beeps are announced three times in succession. This warning is sounded when the meter is switched on.


► Replacing the Battery

To replace the battery, make sure that the meter is turned off.

Step 1: Press the edge of the battery cover and lift it up to remove.

Step 2: Remove the old battery and replace with one 3V CR2032 lithium battery.

Step 3: Close the battery cover. If the battery is inserted correctly, you will hear a “beep” afterwards.

Note:

- Replacing the battery does not affect the test results stored in the memory.
- As with all small batteries, these batteries should be kept away from small children. If swallowed, promptly seek medical assistance.
- Battery might leak chemicals if unused for a long time. Remove the battery if you are not going to use the device for an extended period (i.e., 3 months or more).

- Properly dispose of the battery according to your local environmental regulations.

Caring for Your Device

- To clean the exterior of the device, wipe it with a cloth moistened with tap water or a mild cleaning agent, then dry the device with a soft dry cloth. Do **NOT** rinse with water.
- Do **NOT** use organic solvents to clean the device.

Device Storage

- Storage condition: -20°C to 60°C (-4°F to 140°F), below 95% relative humidity.
- Always store or transport the device in its original storage case.
- Avoid dropping and heavy impact.
- Avoid direct sunlight and high humidity.

Meter Disposal

The used meter should be treated as contaminated and may carry a risk of infection during measurement. The batteries in this used meter should be removed and the meter should be disposed in accordance with local regulations.

The meter falls outside the scope of the European Directive 2002/96/ EC-Directive on waste electrical and electronic equipment (WEEE).

Caring for Your Test Strips

- Storage condition: 4°C to 40°C (39.2°F to 104°F) for glucose test strip, below 85% relative humidity. Do **NOT** freeze.
- Store your test strips in their original vial only. Do not transfer to another container.
- Store test strip packages in a cool and dry place. Keep away from direct sunlight and heat.
- After removing a test strip from the vial, immediately close the vial cap tightly.
- Touch the test strip with clean and dry hands.
- Use each test strip immediately after removing it from the vial.
- Write the opening date on the strip vial label when you first open it. Discard remaining test strips after 6 months.
- Do not use test strips beyond the expiry date. This may cause inaccurate results.
- Do not bend, cut, or alter a test strip in any way.
- Keep the strip vial away from children since the cap and the test strip may be a choking hazard. If swallowed, promptly see a doctor for assistance.

For further information, please refer to the test strip package insert.

Important Control Solution Information

- Use only our control solutions with your device.
- Do not use the control solution beyond the expiry date or 3 months after first opening. Write the opening date on the control solution vial and discard the remaining solution after 3 months.

- It is recommended that the control solution test be done at room temperature 20°C to 25°C (68°F to 77°F). Make sure your control solution, device, and test strips are at this specified temperature range before testing.
- Shake the vial before use, discard the first drop of control solution, and wipe off the dispenser tip to ensure a pure sample and an accurate result.
- Store the control solution tightly closed at temperatures between 2°C to 30°C (35.6°F to 86°F). Do **NOT** freeze.

Reference Values




The device provides you with plasma equivalent blood glucose results.

| Time of day | Normal plasma glucose range for people without diabetes (mg/dL) |
|--------------------------|---|
| Fasting and before meals | < 100 mg/dL (5.6 mmol/L) |
| 2 hours after meals | < 140 mg/dL (7.8 mmol/L) |

Source: *American Diabetes Association (2012). Clinical Practice Recommendations. Diabetes Care, 35 (Supplement 1): S1-100.*

Please consult your doctor to determine a target range that works best for you.

Symbol Information

| Symbol | Referent | Symbol | Referent |
|---|---|---|---|
|  | For in vitro diagnostic use |  | Do not reuse |
|  | Consult instructions for use |  | Temperature limitation |
|  | Use by | | |
|  | Batch code |  | CE Mark |
|  | Serial number |  | Manufacturer |
|  | Keep away from sunlight |  | Dispose of the packaging properly after use |
|  | Keep Dry |  | Caution |
|  | Use within 6 months after opening |  | Sterilized using irradiation |
|  | Disposal of waste equipment |  | Do not use if package is damaged |
|  | 3 Volts DC |  | Humidity limitation |
|  | Model number |  | Battery |
|  | Authorised representative in the European Union |  | Medical Device |

Troubleshooting

If you follow the recommended steps but the problem persists, or error messages other than the ones below appear, please call your local customer service.

Do not attempt to repair the device yourself and never try to disassemble the device under any circumstances.

Result Readings (for glucose test)

| APPEARS | WHEN GLUCOSE | | |
|---|---|---|------------------------------------|
|  | < 20 mg/dL (1.1 mmol/L) | | |
|  | 20-69 mg/dL (1.1-3.8 mmol/L) | | |
|  | AC  | PC  | Gen |
| | 70-129 mg/dL (3.9-7.2 mmol/L) | 70-179 mg/dL (3.9-9.9 mmol/L) | 70-119 mg/dL (3.9-6.6 mmol/L) |
|  | AC  | PC  | Gen |
| | 130-239 mg/dL (7.2-13.3 mmol/L) | 180-239 mg/dL (9.9-13.3 mmol/L) | 120-239 mg/dL (6.7-13.3 mmol/L) |
| KETONE? | ≥ 240 mg/dL (13.3 mmol/L) | | |
|  | > 600 mg/dL (33.3 mmol/L) | | |

Error Message

| Error Message | Cause | What To Do |
|---|--|--|
|  | Appears when the batteries cannot provide enough power for a test. | Replace the batteries immediately. |
|  | Problem with the meter. | Repeat the test with a new test strip. If the meter still does not work, please contact the customer service for assistance. |
|  | Remove the strip after applying blood to the absorbent hole. | Re-test with a new test strip. |
|  | Appears when a used test strip is inserted. | Repeat with a new test strip. |

| | | |
|---|---|--|
|  | Appears when ambient temperature is below system operation range. | System operation range is 10°C to 40°C (50°F to 104°F). Repeat the test after the meter and test strip are in the above temperature range. |
|  | Appears when ambient temperature is above system operation range. | |

Announcement of Error Message by Universal Tone

- Hi or Lo Result: Results larger than 600 mg/dL (33.3 mmol/L) are represented as 999, i.e. three groups of nine short beeps with pauses between the groups.
- Results lower than 20 mg/dL (1.1 mmol/L) are represented as 000, i.e. three long beeps.
- Low Battery Warning: When the batteries are about to be exhausted, 2 quick beeps are announced three times in succession. This warning is sounded when the meter is switched on.
- Other Errors: Other Error Messages are announced by 2 quick beeps four times in succession.

Blood Glucose Measurement

1. If the meter does not display a message after inserting a test strip:

| POSSIBLE CAUSE | WHAT TO DO |
|--|--|
| Batteries exhausted. | Replace the batteries |
| Test strip inserted upside down or incompletely. | Insert the test strip with contact bars end first and facing up. |
| Defective meter or test strips. | Please contact customer services. |

2. If the test does not start after applying the sample:

| POSSIBLE CAUSE | WHAT TO DO |
|---|--|
| Insufficient blood sample. | Repeat the test using a new test strip with larger volume of blood sample. |
| Defective test strip. | Repeat the test with a new test strip. |
| Sample applied after automatic switch-off (2 minutes after last user action). | Repeat the test with a new test strip. Apply sample only when flashing "  " appears on the display. |
| Defective meter. | Please contact customer services. |

3. If the control solution testing result is out of range:

| POSSIBLE CAUSE | WHAT TO DO |
|--|--|
| Error in performing the test. | Read instructions thoroughly and repeat the test again. |
| Control solution vial was poorly shaken. | Shake the control solution vigorously and repeat the test again. |
| Expired or contaminated control solution. | Check the expiry date of the control solution. |
| Control solution that is too warm or too cold. | Control solution, meter, and test strips should be at room temperature (20°C to 25°C / 68°F to 77°F) before testing. |
| Defective test strip. | Repeat the test with a new test strip. |
| Meter malfunction. | Please contact customer services. |

Specifications

| | |
|------------------------------------|---|
| Model No. | FORA COMFORT plus G30a |
| Dimension & Weight | 85(L) x 52(W) x 15(H) mm, 41g |
| Power Source | One CR2032 lithium battery |
| Display | LCD |
| Memory | 450 measurement results with respective date and time |
| External Output | RS232 PC interface |
| Features | Auto electrode insertion detection Auto sample loading detection Auto reaction time count-down Auto switch-off after 3 minutes without action Temperature Warning |
| Operating Conditions | 10°C to 40°C (50°F to 104°F), below 85% R.H. (non-condensing) |
| Storage / Transportation Condition | -20°C to 60°C (-4°C to 140°F), below 95% R.H. |
| Measurement Units | mg/dL or mmol/L |
| Measurement Range | 20 to 600mg/dL (1.1 to 33.3mmol/L) |

This device has been tested to meet the electrical and safety requirements of: IEC/EN 61010-1, IEC/EN 61010-2-101, EN 61326-1, IEC/EN 61326-2-6, EN 301 489-17, EN 300 328.

WARRANTY TERMS AND CONDITIONS

With respect to disposable products, ForaCare Suisse warrants to the original purchaser that, at time of delivery, each standard product manufactured by ForaCare Suisse shall be free from defects in material and workmanship and, when used for the purposes and indications described on the labeling, is fit for the purposes and indications described on the labeling. All warranties for a product shall expire as of the product expiration date, or if none, after two (2) year from the original date of purchase, as long as it has not been modified, altered, or misused. ForaCare Suisse warranty hereunder shall not apply if:

(i) a product is not used in accordance with its instructions or if it is used for a purpose not indicated on the labeling; (ii) any repairs, alterations or other work has been performed by the buyer or others on such item, other than work performed with ForaCare Suisse's authorisation and according to its approved procedures; or (iii) the alleged defect is a result of abuse, misuse, improper maintenance, accident or the negligence of any party other than ForaCare Suisse. The warranty set forth herein is conditioned upon proper storage, installation, use and maintenance in accordance with applicable written recommendations from ForaCare Suisse.

The warranty furnished hereunder does not extend to damaged items purchased hereunder resulting in whole or in part from the use of components, accessories, parts or

supplies not furnished by ForaCare Suisse.

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ForaCare Suisse AG
MedNet EC-REP GmbH