

FORA[®] ADVANCED pro GD40

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














(f)

Safety Information

Read the following *Safety Information* thoroughly before using the device.

- Use this device **ONLY** for the intended use described in this manual.
- Do **NOT** use accessories which are not specified by the manufacturer.
- Do **NOT** use the device if it is not working properly or damaged.
- This device does **NOT** serve as a cure for any symptoms or diseases. The data measured is for reference only. Always consult your doctor to have the results interpreted.
- The blood glucose test strip can be used for the testing of newborns; The β -Ketone test strip must **NOT** be used for the testing of newborns.
- Before using this device to test blood glucose or β -Ketone, read all instructions thoroughly and practice the test. Carry out all the quality control checks as directed.
- Keep the device and testing supplies away from young children. Small items such as the battery cover, batteries, test strips, lancets and vial caps are choking hazards.
- Use of this instrument in a dry environment, especially if synthetic materials are present (synthetic clothing, carpets etc.) may cause damaging static discharges that may cause erroneous results.
- Do **NOT** use this instrument in close proximity to sources of strong electromagnetic radiation, as these may interfere with the correct operation.
- Proper maintenance and timely calibration with the control solution are essential to the longevity of your device. If you are concerned about the accuracy of the measurement, please contact the place of purchase or customer service representative for assistance.

KEEP THESE INSTRUCTIONS IN A SAFE PLACE

Important Information

- Severe dehydration and excessive water loss may cause readings which are lower than actual values. If you believe you are suffering from severe dehydration, consult a healthcare professional immediately.
- If your blood glucose or β -Ketone results are lower or higher than usual, and you do not have symptoms of illness, first repeat the test. If you have symptoms or continue to get results higher or lower than usual, follow the treatment advice of your healthcare professional.
- Use only fresh whole blood sample to test your blood glucose or β -Ketone. Using other substances will lead to incorrect results.
- If you are experiencing symptoms that are inconsistent with your blood glucose or β -Ketone test results and you have followed all instructions described in this owner's manual, call your healthcare professional.
- We do not recommend using this product on severely hypotensive individuals or patients in shock. Readings which are lower than actual values may occur for individuals experiencing a hyperglycaemic-hyperosmolar state, with or without ketosis. Please consult the healthcare professional before use.
- The measurement unit used for indicating the concentration of blood or plasma glucose can have mg/dL or mmol/L. The approximate calculation rule for conversion of mg/dL in mmol/L is:

mg/dL	Divided by 18	= mmol/L
mmol/L	Times 18	= mg/dL

For example:

1) $120 \text{ mg/dL} \div 18 = 6.6 \text{ mmol/L}$

2) $7.2 \text{ mmol/L} \times 18 = 129 \text{ mg/dL}$ approximately

- The HCT is the volume percentage (%) of red blood cells in the blood sample. The Hb is the haemoglobin (g/dL) levels in the blood sample. The approximate calculation rule for Hb is: $\text{Hb} = \text{HCT} \times 0.34$ (Note: Haematocrit =HCT; Haemoglobin=Hb).
- The measurement unit used for indicating the concentration of β -Ketone is mmol/L.

Introduction

Intended Use

This system is intended for use outside the body (in vitro diagnostic use) by people with diabetes at home and by health care professionals in clinical settings as an aid to monitor the effectiveness of diabetes control. It is intended to be used for the quantitative measurement of glucose (sugar) with capillary from finger tip, venous, arterial, and to measure β -hydroxybutyrate (Ketone) in fresh whole blood samples from the finger. It should not be used for the diagnosis of diabetes, or screening for diabetes mellitus.

The test for hematocrit (HCT) as part of the system, is intended for use in the in vitro quantification of packed red blood cell volume fraction in capillary whole blood as an aid in monitoring the status of total volume of red blood cells. The reading of HCT is only for the determination of whether the blood test sample is within the acceptable range of Blood Glucose Plus β -Ketone monitoring system. It should not be used for the diagnosis of anemia or erythrocytosis.

Professionals may test blood glucose and hematocrit with capillary sampling from a finger tip or with venous, arterial, and neonatal blood from the heel. Test β -Ketone with a capillary sampling from the finger tip. Use only heparin for anticoagulation of whole blood.

Home use is limited to capillary blood from the finger tip.

Test Principle

Your system measures the amount of glucose or β -Ketone in whole blood. The glucose or β -Ketone testing is based on the measurement of electrical current generated by the reaction of glucose or β -Ketone with the reagent of the strip. The meter measures the current, calculates the blood glucose or β -Ketone level, and displays the result. The strength of the current produced by the reaction depends on the amount of glucose or β -Ketone in the blood sample.

Product Overview (a)

- | | |
|-----------------------|------------------------------------|
| 1. Test Slot | 5. Data Port (for GD40g) |
| 2. MAIN Button | Bluetooth Switch (for GD40h) |
| 3. Display Screen | 6. Bluetooth Indicator (for GD40h) |
| 4. Test Strip Ejector | 7. SET Button |
| | 8. Battery Compartment |

Screen Display (b)

- | | |
|-----------------------------------|---------------------------|
| 1. Blood Drop Symbol | 9. Measurement Mode |
| 2. Code | 10. Measurement Unit |
| 3. Ketone Warning / Ketone Symbol | 11. HCT Symbol |
| 4. Test Result | 12. HCT Result |
| 5. Warning Symbol | 13. Hb Symbol |
| 6. Low Battery Symbol | 14. Date/Time/Hb Result |
| 7. Day Average | 15. Test Strip Symbol |
| 8. Memory Symbol | 16. Universal Tone Symbol |

Getting Started

Initial Setup


Before using the device for the first time or after you have changed the batteries, enter the following initial setup process:

Step 1: Enter the Setting Mode

Open the battery cover and press **SET**. The display turns on.

Step 2: Configuring the Settings (Date, Time Format, Time, Universal Tone, Memory Deletion)

Press **MAIN** repeatedly to adjust the value or enable/disable the setting. Then press **SET** to confirm the setting and switch to another field.

Note:

- You may set up Beep On, Universal Tone On and Beep Off by pressing M to change it. (v)
- When Universal Tone is turned on, the meter guides you through the blood glucose test in a series of beep tones; it also outputs the result as a series of beeps.
- During memory deletion, select 'No' to keep all saved results.
- If the device is idle for 3 minutes during the setting mode, it will turn off automatically.

Before Testing

Calibration

You must calibrate the device every time you begin to use a new box of β -Ketone test strips by setting the meter with the correct code. To ensure test accuracy, make sure the code number displayed on the display screen matches the number printed on the test strip vial or individual foil pack.

How to Code Your Meter

1. Insert the code strip into the test slot of the device. Wait for the device to display the code number.

Important!

Make sure the code number on display, code strip, and test strip vial or individual foil pack are the same. **(c)**

2. Remove the code strip, the display will show "OFF" indicating the device has finished coding and is ready for β -Ketone or blood glucose testing.

Checking the Code Number

Make sure that the code number displayed on the device matches the number on the test strip vial or individual foil pack before you proceed. **(d)**

If it matches, you can proceed with your test. Otherwise, please stop testing and repeat the calibration procedure.

If the problem persists, contact Customer Service for further assistance.

Important!

• It is important to make sure that the LCD displayed code is the same as the code on the test strip vial or individual foil pack before testing. Failure to do so will cause inaccurate results.

Control Solution Testing

Our Control Solution contains a known amount of glucose or β -Ketone that reacts with test strips and is used to ensure your device and test strips are working together correctly.

Do a control solution test when:

- ✓ you first receive the device.
- ✓ at least once a week to routinely check the device and test strips.
- ✓ you begin using a new vial of test strips.
- ✓ you suspect the device or test strips are not working properly.
- ✓ your blood glucose or β -Ketone test results are not consistent with how you feel, or if you think the results are not accurate.
- ✓ you have dropped or think you may have damaged the device.

To perform the control solution test, do the following:

1. Insert the test strip into the test slot of the device. Wait for the device to display the test strip and blood drop symbol. **(e)**
2. The meter will detect the difference between control solution and blood samples automatically. It will automatically mark the result as a control solution test with "QC" display.
3. Apply the control solution.

Shake the control solution vial thoroughly before use. Squeeze out a drop and wipe it off, then squeeze another drop and place it on the tip of the vial cap. Hold the device to move

the absorbent hole of the test strip to touch the drop. Once the confirmation window is filled completely, the device will begin counting down. **(g)**

Note:

- To avoid contaminating the control solution, do not directly apply the control solution onto a strip.

4. Read and compare the result. (h)

After counting down to 0, the test result of the control solution will appear on the display. Compare this result with the range printed on the test strip vial or individual foil pack and it should fall within this range. If the test result is out of range, read the instructions again and repeat the control solution test.

Note:


- Control solution test results are stored in the memory.
- The control solution range printed on the test strip vial or individual foil pack is for control solution use only. It is not a recommended range for your blood glucose level.
- Refer to the **Maintenance** section for important information about your control solutions.

Out-of-range results:

- If you continue to have test results falling outside the range printed on the test strip vial, the meter and strips may not be working properly.
- There is no HCT and HB display functions while the meter is in QC mode.

Testing Your Blood Glucose


Test Strip Appearance

1. Blood Glucose Test Strip

2. Blood Glucose / Hematocrit / Hemoglobin Test Strip

3. Blood β -Ketone Test Strip

Inserting a Test Strip

Insert the test strip into its slot.

Important!

The front side of test strip should face up when inserting test strip. Test results might be wrong if the contact bar is not fully inserted into the test slot.

Important!

To reduce the chance of infection:

- Never share a lancet or a lancing device.
- Always use a new, sterile lancet. Lancets are for single use only.
- Avoid getting hand lotion, oils, dirt, or debris in or on the lancets and the lancing device.

Preparing the Lancing Device

1. Remove the cap. **(i)**
2. Insert a new lancet firmly into the white lancet holder cup.
3. Remove the protective disk on the lancet. **(j)**
Hold the lancet firmly in place and twist off the protective disk.
4. Replace the cap until it snaps or clicks into place.
5. Rotate the dial to set the desired lancing depth. **(k)**
6. Pull the cocking control out until the orange bar appears on the release button window. **(l)**

Obtaining the Blood Sample

Please follow the suggestions below before obtaining a drop of blood:

- Wash and dry your hands before starting.
- Select the puncture site at fingertips.
- Clean the puncture site using cotton moistened with 70% alcohol and let it air dry.
- Rub the puncture site for about 20 seconds before penetration.

❖ Blood from the fingertip

1. Press the lancing device tip firmly against the lower side of your fingertip.
2. Press the release button to prick your finger. A click indicates that the puncture is complete. **(m)**

Note:

- Choose a different spot each time you test. Repeated punctures at the same spot may cause soreness and calluses.
- It is recommended to discard the first drop of blood as it might contain tissue fluid, which may affect the test result.

Make Blood Glucose or β -Ketone Test

1. Insert the test strip into its slot.
Wait for the device to display the test strip "", blood drop "", and "**GLU**" or "**KETONE**".

2. Select the appropriate measuring mode by pressing **MAIN**.

- General Tests (**Gen**) - any time of day without regard to time since last meal.
- AC (**AC**) - no food intake for at least 8 hours.
- PC (**PC**) - 2 hours after a meal.

3. Obtain a blood sample. (**o**)

Use the pre-set lancing device to puncture your desired site. After penetration, discard the first drop of blood with a clean tissue or cotton. Gently squeeze the punctured area to obtain another drop of blood. Be careful NOT to smear the blood sample. The volume of blood sample must be at least 1.1 microliter (μL) for blood glucose or at least 1.0 microliter (μL) for β -Ketone.

4. Apply the blood sample. (**p**)

Hold the blood drop to touch the absorbent hole of the test strip. Blood will be drawn in and after the confirmation window is completely filled, the device begins counting down. Do not remove your finger until you hear a beep sound.

5. Read your result. (**q**)

The results of your blood glucose test with HCT/ Hb levels or the result of your β -Ketone test will appear after the meter counts down to 0. The results will be stored automatically in the meter memory.

Announcement of Test Result by Universal Tone

The blood glucose results will be broken down into individual digits and each digit represents the corresponding number of beeps.

The result is announced three times in succession and each time is preceded by two quick beeps. So you will hear: 2 quick beeps – results – 2 quick beeps – results – 2 quick beeps – result.

For mg/dL meters, the hundreds are always announced, even when the result is below 100.

Examples:

80 mg/dL is announced as 1 long beep (0) – 1 single pause – 8 single beeps (8) – 1 single pause – 1 long beep (0)

182 mg/dL is announced as 1 single beep (1) – 1 single pause – 8 single beeps (8) – 1 single pause – 2 single beeps (2)

For mmol/L meters, the tens are always announced, even when the result is below 10. The decimal point is represented by 1 quick beep.

Examples:

6.0 mmol/L is announced as 1 long beep (0) – 1 single pause – 6 single beeps (6) – 1 single pause – 1 quick beep (.) – 1 single pause – 1 long beep (0)

Note:

Information or warnings in the form of symbols displayed together with the results are not announced acoustically.

Disposing Used Test Strip and Lancet

To remove the used test strip, simply press the **Test Strip Ejector** to eject the used test strip. The device will automatically turn off after the test strip is ejected.

To remove the used lancet, remove the lancet from the lancing device after you have finished testing.

Discard your used strip and lancet properly in a puncture resistant container.

Important!

The used lancet and test strip may be biohazards. Please consult your health care provider for proper disposal complying with your local regulations.

Reviewing Test Results

Your device stores the 1000 most recent test results along with the respective dates and times in its memory.

To review all test results, do the following:

1. Press and release **MAIN**. The “” icon appears on the screen.
2. Press **MAIN** to review the test results stored in the device. Press **MAIN** repeatedly to review other test results stored in the device. **(r)** After the last test result, press **MAIN** again and the device will be turned off.

To review the day-average test results, do the following:

1. Press and hold **MAIN** for 3 seconds until the “**DAY AVG**” icon appears. Release **MAIN** and then your 7-day average result

measured in general mode will appear on the display.

2. Press **MAIN** to review 14, 21 and 30 day average results stored in each measuring mode in the order of Gen, AC, then PC. **(s)**

Note:

- Press and hold **MAIN** for 5 seconds to exit the memory mode or leave it without any action for 3 minutes. The device will turn off automatically.
- If using the device for the first time, the “---” icon will appear when you recall the test results or review the average result. This indicates that there is no test result in the memory.
- Control solution results are **NOT** included in the day average.

Announcement of Memory Result by Universal Tone

Only the most recent result that was saved can be announced acoustically. If you press the M button to turn the meter on, you will first hear the Long-Beep which stands for power-on and then the most recent result.

Only the average for the last 7 days is announced acoustically. If the 7-day average cannot be calculated, three horizontal bars are displayed. This is signaled acoustically with 3 long beeps representing 3 zeroes.

Appendix

Transferring Data

There are two types of connection (either USB/RS-232 or Bluetooth) to transfer results. Please check your device box for the transmission method of your device.

1. Data Transmission Via USB/RS-232 Cable (For GD40g only)

i. Install the software on your computer

Download Health Care System Software and instruction manual provided on the ForaCare Suisse AG Web site: <http://www.foracare.ch>. Follow the instructions to install the software on your computer.

ii. Connect the device with your computer using a USB/RS-232 cable

“PC” will appear on the display, indicating that the device is ready to transmit data. **(t)**

iii. Transfer data to your computer

Follow the software on-screen instructions to transmit data. The data transmitted will include results with date and time. Remove the cable and the device will automatically turn off.

2. Data Transmission Via Bluetooth (For GD40h only)

i. Install the software on your computer

Download Health Care System Software and instruction manual provided on the ForaCare Suisse AG Web site:

<http://www.foracare.ch>. Follow the instructions to install the software on your computer.

ii. Connect the device with your computer

With the device turned off, slide the **Bluetooth switch** on the device to initiate the data transmission. The Bluetooth indicator flashes in blue.

iii. Transfer data to your computer

Follow the software on-screen instructions to transmit data. The data transmitted will include results with date and time.

Note:

- Make sure your PC has Bluetooth capability before transmitting the data and the device is within the receiving range.
- When the device is connected to the PC, it is unable to perform any blood glucose or β -Ketone test.

Maintenance

Changing Battery

When the battery is low, the “” icon appears with E-b. **(u)**

To change the batteries, do the following:

1. Press the edge of the battery cover and lift it up to remove the cover.
2. Remove the old batteries and replace them with two 1.5V AA size alkaline batteries.

3. Close the battery cover.

When you hear a “beep” sound, it indicates that the batteries are inserted correctly.

CAUTION
RISK OF EXPLOSION IF BATTERY IS REPLACED BY AN
INCORRECT TYPE.
DISPOSE OF USED BATTERIES ACCORDING TO THE
INSTRUCTIONS.

Note:

- Replacing the batteries does not affect the test results stored in the memory.
- Keep these batteries away from small children. If swallowed, promptly seek medical assistance.
- Batteries may leak chemicals if unused for a long time. Remove the batteries if you are not going to use the device for an extended period of time.
- Properly dispose of the used batteries according to your local environmental regulations.

Caring for Your Device

- To clean the exterior of the device, wipe it with a cloth moistened with tap water or a mild cleaning agent, then dry the device with a soft dry cloth. Do **NOT** rinse with water.
- Do **NOT** use organic solvents to clean the device.

Device Storage

- Storage condition: -20°C to 60°C (-4°F to 140°F), below 95% relative humidity.
- Always store or transport the device in its original storage case.

- Avoid dropping and heavy impact.
- Avoid direct sunlight and high humidity.

Meter Disposal

The used meter should be treated as contaminated and may carry a risk of infection during measurement. The batteries in this used meter should be removed and the meter should be disposed in accordance with local regulations.

The meter falls outside the scope of the European Directive 2002/96/ EC-Directive on waste electrical and electronic equipment (WEEE).

Caring for Your Test Strips

- Storage condition: 2°C to 32°C (35.6°F to 89.6°F) for glucose test strip and 4°C to 30°C (39.2°F to 86°F) for β -Ketone test strip, below 85% relative humidity. Do **NOT** freeze.
- Store your test strips in their original vial only. Do not transfer to another container.
- Store test strip packages in a cool and dry place. Keep away from direct sunlight and heat.
- After removing a test strip from the vial, immediately close the vial cap tightly.
- Touch the test strip with clean and dry hands.
- Use each test strip immediately after removing it from the vial.
- Write the opening date on the strip vial label when you first open it. Discard remaining test strips after 6 months.
- Do not use test strips beyond the expiry date. This may cause inaccurate results.
- Do not bend, cut, or alter a test strip in any way.

- Keep the strip vial away from children since the cap and the test strip may be a choking hazard. If swallowed, promptly see a doctor for help.

For further information, please refer to the test strip package insert.

Important Control Solution Information

- Use only our control solutions with your device.
- Do not use the control solution beyond the expiry date or 3 months after first opening. Write the opening date on the control solution vial and discard the remaining solution after 3 months.
- It is recommended that the control solution test be done at room temperature 20°C to 25°C (68°F to 77°F). Make sure your control solution, device, and test strips are at this specified temperature range before testing.
- Shake the vial before use, discard the first drop of control solution, and wipe off the dispenser tip to ensure a pure sample and an accurate result.
- Store the control solution tightly closed at temperatures between 2°C to 30°C (35.6°F to 86°F). Do **NOT** freeze.

Reference Values

- The device provides you with plasma equivalent blood glucose results.

Time of day	Normal plasma glucose range for people without diabetes (mg/dL)
Fasting and before meal	< 100 mg/dL (5.6 mmol/L)
2 hours after meals	< 140 mg/dL (7.8 mmol/L)

Source: American Diabetes Association (2012). *Clinical Practice Recommendations. Diabetes Care, 35 (Supplement 1): S1-100.*

- Normal Hematocrit test results may vary, but in general they are as follows:

- Male: 40.7 to 50.3%
- Female: 36.1 to 44.3%

If the hematocrit test result is more than 70%, or hemoglobin is more than 23.8 (g/dL), the blood glucose test may be invalid. Please redo the test, and consult your doctor if the result is repeatedly invalid.

- The β -Ketone test measures Beta-Hydroxybutyrate (β -OHB), the most important of the three β -Ketone bodies in the blood. Normally, levels of β -OHB are expected to be less than 0.6 mmol/L. β -OHB levels may increase if a person fasts, exercises vigorously or has diabetes and becomes ill. If your β -Ketone result is 0.0 mmol/L, repeat the β -Ketone test with new test strips. If the same message appears again or the result does not reflect how you feel, contact your healthcare professional. Follow your healthcare professional's advice before you make any changes to your diabetes medication programme. If your β -Ketone result is between 0.6 and 1.5 mmol/L, this may indicate a development of a problem that may require medical assistance. Follow your healthcare professional's instructions. If your β -Ketone result is higher than 1.5 mmol/L, contact your healthcare professional promptly for advice and assistance. You may be at risk of developing diabetic ketoacidosis (DKA).

Please consult your doctor to determine a target range that works best for you.

Symbol Information

Symbol	Referent	Symbol	Referent
	For in vitro diagnostic use		Do not reuse
	Consult instructions for use		Temperature limitation
	Use by		CE Mark
	Batch code		Manufacturer
	Serial number		Dispose of the packaging properly after use
	Keep away from sunlight		Caution
	Keep Dry		Sterilized using irradiation
	Use within 6 months after opening		Do not use if package is damaged
	Humidity limitation		Authorised representative in the European Community
	3 Volts DC		Disposal of waste equipment
	Medical device		Model No.

Troubleshooting

If you follow the recommended action but the problem persists, or error messages other than the ones below appear, please call your local customer service.

Result Readings (for glucose test)

Message	What it Means
 <p>Gen LO HCT 58% Hb 19.72g/dL</p>	< 10 mg/dL (0.5 mmol/L)
 <p>Gen KETONE 300 HCT 58% Hb 19.72g/dL</p>	≥ 240 mg/dL (13.3 mmol/L)
 <p>Gen HI HCT 58% Hb 19.72g/dL</p>	> 600 mg/dL (33.3mmol/L)

Result Reading (for β -Ketone test)

Message	What it Means
Lo	< 0.1 mmol/L
06	0.1 to 8.0 mmol/L
Hi	> 8.0 mmol/L

Error Message

Error Message	Cause	What To Do
E-b	The batteries cannot provide enough power for a test.	Replace the batteries immediately.
E-U	Strip has been used.	Repeat the test with a new strip.
E-E E-3	Problem in operation.	Review the instructions and repeat the test with a new strip. If problem persists, contact the local customer service for assistance.
E-0 E-A		
E-F	You may have removed the strip after applying blood, or insufficient blood volume.	Review the instructions and repeat test with a new test strip.
	HCT>75%	

E-C	Wrong code strip is inserted or other coding errors.	Make sure to insert the right code strip for the meter. Check if the code numbers on the code strip and the strip vial label are the same.
E-2	Expired code strip.	Repeat the test with a new code strip.
E-t	Ambient temperature is below the system's operation range.	System operational range is 10°C to 40°C (50°F to 104°F). Repeat the test after the device and test strip have reached the above temperature.
E-t	Ambient temperature is above the system's operation range.	

Announcement of Error Message by Universal Tone

- Hi or Lo Result: Results larger than 600 mg/dL (33.3 mmol/L) are represented as 999, i.e. three groups of nine short beeps with pauses between the groups.
- Results lower than 10 mg/dL (0.5 mmol/L) are represented as 000, i.e. three long beeps.
- Low Battery Warning: When the batteries are about to be exhausted, 2 quick beeps are announced three times in succession. This warning is sounded when the meter is switched on.
- Other Errors: Other Error Messages are announced by 2 quick beeps four times in succession.

Blood Glucose Measurement

Symptom	Cause	What To Do
The device does not display a message after inserting a test strip.	Batteries exhausted.	Replace the batteries.
	Test strip inserted upside down or incompletely.	Insert the test strip with contact bars end first and facing up.
	Defective device or test strips.	Please contact customer service.

The test does not start after applying the sample.	Insufficient blood sample.	Repeat the test using a new test strip with larger volume of blood sample.
	Defective test strip.	Repeat the test with a new test strip.
	Sample applied after the device is automatically turned off.	Repeat the test with a new test strip. Apply sample only when flashing "⚡" appears on the display.
	Defective device.	Please contact customer service.
The control solution testing result is out of range.	Error in performing the test.	Read instructions thoroughly and repeat the test again.
	Control solution vial was poorly shaken.	Shake the control solution vigorously and repeat the test again.
	Expired or contaminated control solution.	Check the expiration date of the control solution.
	Control solution that is too warm or too cold.	Control solution, device, and test strips should be at room temperature (20°C to 25°C / 68°F to 77°F) before testing.
	Defective test strip.	Repeat the test with a new test strip.
	Device malfunction.	Please contact customer service.

Specifications

Memory	1000 measurement results with respective date and time
Dimensions	110 (L) x 57 (W) x 25 (H) mm
Power Source	3V  (Two 1.5V AAA alkaline batteries)
Weight	71g
External output	Serial Cable (for GD40g) Bluetooth (for GD40h)
Features	Auto electrode insertion detection Auto sample loading detection Auto reaction time count-down Auto switch-off after 3 minutes without action Temperature warning
Operating Condition	10°C to 40°C (50°F to 104°F), below 85% R.H. (non-condensing)
Storage/Transportation Condition	-20°C to 60°C (-4°F to 140°F), below 95% R.H
Measurement Units	either mg/dL or mmol/L for glucose test and fixed mmol/L for β -Ketone test
Measurement Range	10 to 600 mg/dL (0.5 to 33.3 mmol/L) for glucose test and 0.1~8.0 mmol/L for β -Ketone test
Hematocrit range	0~70% for glucose testing and 10~70% for β -Ketone test
Test Sample	Capillary whole blood
Test Result	Glucose measurements are reported as plasma equivalents

This device has been tested to meet the electrical and safety requirements of: IEC/EN 61010-1, IEC/EN 61010-2-101, EN 61326-1, IEC/EN 61326-2-6, EN 301 489-17, EN 300 328.

WARRANTY TERMS AND CONDITIONS

With respect to disposable products, ForaCare Suisse warrants to the original purchaser that, at time of delivery, each standard product manufactured by ForaCare Suisse shall be free of defects in material and workmanship and, when used for the purposes and indications described on the labeling, is fit for the purposes and indications described on the labeling. All warranties for a product shall expire as of the product expiration date, or if none, after two (2) year from the original date of purchase, as long as it has not been modified, altered, or misused. ForaCare Suisse warranty hereunder shall not apply if:

(i) a product is not used in accordance with its instructions or if it is used for a purpose not indicated on the labeling; (ii) any repairs, alterations or other work has been performed by the buyer or others on such item, other than work performed with ForaCare Suisse's authorisation and according to its approved procedures; or (iii) the alleged defect is a result of abuse, misuse, improper maintenance, accident or the negligence of any party other than ForaCare Suisse. The warranty set forth herein is conditioned upon proper storage, installation, use and maintenance in accordance with applicable written recommendations of ForaCare Suisse.

The warranty furnished hereunder does not extend to damaged items purchased hereunder resulting in whole or in part from the use of components, accessories, parts or supplies not furnished by ForaCare Suisse.

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MedNet EC-REP GmbH